


C O L O U R S		RAZOR BLADE		ORDER FORM		Ver. 9.0 Rev-072513			
DEALER				COLOURS INFORMATION For Office Use Only					
Dealer: Discover Your Mobility				Order Date:		S/N:			
Account:		P.O.#:		RTS					
Ordered By:				RTS:		Phone:			
Address: 66 East 10 Mile Rd				Email:					
City: Hazel Park		State: MI		Zip: 48030		SHIP TO			
Phone: (866) 868-9694		Fax:				Name:		Ship Via:	
Email:				Address:					
Customer Name:				City:		State:		Zip:	
Customer Weight:				Phone:		Fax:			
CHAIR INFORMATION					PRICE				
RAZOR BLADE K0005 Fixed FSH-RSH Max. Capacity 250Lbs NOTE <u>Main Frame : ZK55 & ZK60 High Grade A7000 Series Aluminum</u> <u>Main Frame : Lifetime Warranty</u> *See warranty policy					Base Wheel Chair Price K0005 \$ 1,920				
					Total Section A \$				
					Total Section B \$				
					Total Section C \$				
					Total Section D \$				
					Total Section E \$				
					Total Section F (Custom Parts Color) \$				
					ORDER GRAND TOTAL \$				

A FRAME COLORS & PARTS COLORS


FRAME COLORS			
Powder-Co	☐ Gloss Black	☐ Twilight Flake	☐ Yellow
	☐ Candy Blue	☐ Candy Purple	☐ Candy Green
	☐ Matte Gray	☐ Matte Purple	☐ Matte Blue
	☐ Neon Red	☐ Neon Green	☐ Orange Sparkle
		☐ Granny Sparkle Green	☐ Sparkle Sky Blue
		☐ Vermilion	☐ Reed Green
		☐ Blue Lilac	STD
Special	☐ Polished		
	☐ 24 karat Gold Plated	\$ 800	
PARTS COLORS (Custom Parts Color Section, Please see Section F)			
STD Color	Black Color (*Suspension Springs-Silver Color		\$ STD
Custom Color Main Bolt on Parts	One Color Package		
Including: Backrest, Foot Plate, Footrest Rings, Backrest Brackets, Saddle Stops,	☐ Powder-Coated ()	\$ 500	
Forks & Knurl Nuts, Caps, Front Suspension Springs, Axle Tube, Axle Clamps, Axle Plates,	☐ Polished	\$ 750	
Axle Bases, Camber Discs ☆marked parts	☐ 24 karat Gold Plated	\$ CFP	
Total Section A			\$

B FRAME MEASUREMENTS 1/2



1. SEAT WIDTH	☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" 	\$	STD								
2. SEAT DEPTH * Only Choose One Option	☐ Net Seat Depth ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19" 	\$	STD								
	☐ 20" E2203	\$	300								
	☐ Depth Adj. Frame ☐ Short Frame (16" Depth) Set at ☐ 13" ☐ 14" ☐ 15" ☐ 16"	\$	STD								
	☐ Long Frame (19" Depth) Set at ☐ 16" ☐ 17" ☐ 18" ☐ 19"	\$	STD								
	☐ Positioning Frame Add (")	\$	150								
	☐ Ergo Seat Tubing Depth (") *Measurement must be 4" longer than your Center of Gravity	\$	200								
3. FOOTREST WIDTH	☐ STD (Seat Width Minus 5") ☐ Taper on Top ☐ Taper on Bottom	\$	STD								
	☐ Custom (") ☐ Taper on Top ☐ Taper on Bottom	\$	150								
	☐ No Taper (Seat Width Minus 2")	\$	STD								
FOOT PLATE STYLE	☐ Standard single loop with adjust. height & angle. STD Foot Plate K0040	\$	STD								
	☐ Single Flip-Back-single loop (no riser avail) with adjust. height & angle. STD Foot Plate K0040	\$	150								
	STD Foot Plate Size Chart <table border="1" style="display: inline-table; vertical-align: top;"> <thead> <tr> <th>Your Footrest Width</th> <th>Foot Plate Size</th> </tr> </thead> <tbody> <tr> <td>8.5"-10.5"</td> <td>W6.5"xD5"</td> </tr> <tr> <td>11.0"-13.0"</td> <td>W8.5"xD6.5"</td> </tr> <tr> <td>13.5"-17.0"</td> <td>W11"xD6.5"</td> </tr> </tbody> </table>	Your Footrest Width	Foot Plate Size	8.5"-10.5"	W6.5"xD5"	11.0"-13.0"	W8.5"xD6.5"	13.5"-17.0"	W11"xD6.5"		
Your Footrest Width	Foot Plate Size										
8.5"-10.5"	W6.5"xD5"										
11.0"-13.0"	W8.5"xD6.5"										
13.5"-17.0"	W11"xD6.5"										
CUSTOM FOOT PLATE	☐ Custom Foot Plate. with adjust. height & angle K0040 *Available Solid Plate Only ☐ Drill 2Sets Extra Holes in Footplate	\$	75								
	☐ Custom Tubular Footrest	\$	250								
	☐ With ABS Foot Plate (\$65 up charge)										
4. FRAME ANGLE	☐ 75° ☐ 80° ☐ 85° ☐ 90° to the ground	\$	STD								
	☐ Custom (°) to the ground	\$	NCO								

B FRAME MEASUREMENTS 2/2

5. BACKREST HEIGHT * Rigidizer bar 3.25" deep std.		☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18"	\$ STD														
TYPE	☐ Locking mechan. with adjust. height & angle		\$ STD														
	☐ Double Locking mechan. with adjust. height & angle		\$ 25														
	☐ Folding with adjust. height & angle		\$ NCO														
	☐ Welded Backrest Angle Set at (") to the ground *Non Folding Backrest only ..		\$ 150														
CUSTOM	☐ Welded Side Guards *Non Folding Backrest only		\$ 450														
	☐ Custom Backrest (ie: double rigidizer bars, custom depth, custom height, etc) 		\$ 150														
6. BACKREST ANGLE		Set at (°) to the ground *STD 90°	\$ STD														
7. FRONT SEAT HEIGHT		☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19" ☐ 19.5" ☐ 20" ☐ 20.5" ☐ 21"	\$ STD														
STD Front Seat Height Ranges <table border="1" data-bbox="532 800 1005 913"> <thead> <tr> <th>Caster Wheel Size</th> <th>Front Seat Height</th> </tr> </thead> <tbody> <tr> <td>3" Caster</td> <td>16.0" - 19.5"</td> </tr> <tr> <td>4" Caster</td> <td>16.0" - 20.0"</td> </tr> <tr> <td>5" Caster</td> <td>16.0" - 21.0"</td> </tr> </tbody> </table>				Caster Wheel Size	Front Seat Height	3" Caster	16.0" - 19.5"	4" Caster	16.0" - 20.0"	5" Caster	16.0" - 21.0"						
Caster Wheel Size	Front Seat Height																
3" Caster	16.0" - 19.5"																
4" Caster	16.0" - 20.0"																
5" Caster	16.0" - 21.0"																
☐ CUSTOM FRONT SEAT HEIGHT (")		\$ 250															
8. REAR SEAT HEIGHT		☐ 14.5" ☐ 15.0" ☐ 15.5" ☐ 16.0" ☐ 16.5" ☐ 17.0" ☐ 17.5" ☐ 18.0" ☐ 18.5" ☐ 19.0" 	\$ STD														
STD Rear Seat Height Ranges *No Custom Modifications <table border="1" data-bbox="532 1031 1005 1228"> <thead> <tr> <th>Rear Tire Size</th> <th>Rear Seat Height</th> </tr> </thead> <tbody> <tr> <td>24"x1"(540)</td> <td>14.5" - 18.0"</td> </tr> <tr> <td>24"x1 3/8"(540)</td> <td>15.0" - 18.0"</td> </tr> <tr> <td>25"x1"(559)</td> <td>15.5" - 18.0"</td> </tr> <tr> <td>24"x2"(540)</td> <td>15.5" - 18.5"</td> </tr> <tr> <td>26"x1"(590)</td> <td>15.5" - 19.0"</td> </tr> <tr> <td>25"x1.95"(559)</td> <td>16.0" - 19.0"</td> </tr> </tbody> </table>				Rear Tire Size	Rear Seat Height	24"x1"(540)	14.5" - 18.0"	24"x1 3/8"(540)	15.0" - 18.0"	25"x1"(559)	15.5" - 18.0"	24"x2"(540)	15.5" - 18.5"	26"x1"(590)	15.5" - 19.0"	25"x1.95"(559)	16.0" - 19.0"
Rear Tire Size	Rear Seat Height																
24"x1"(540)	14.5" - 18.0"																
24"x1 3/8"(540)	15.0" - 18.0"																
25"x1"(559)	15.5" - 18.0"																
24"x2"(540)	15.5" - 18.5"																
26"x1"(590)	15.5" - 19.0"																
25"x1.95"(559)	16.0" - 19.0"																
9. CENTER OF GRAVITY		☐ 1.0" ☐ 1.5" ☐ 2.0" ☐ 2.5" ☐ 3.0" ☐ 3.5" ☐ 4.0" ☐ 4.5" ☐ 5.5" ☐ 6.0" ☐ 6.5" ☐ 7.0"	\$ STD														
10. CAMBER ANGLE		☐ 2° ☐ 4° ☐ 6° 	\$ STD														
11. FRONT OF SEAT SLING TO FRONT OF FOOTREST		★ ☐ 15.0" ★ ☐ 15.5" ★ ☐ 16.0" ☐ 16.5" ☐ 17.0" ☐ 17.5" ☐ 18.0" ☐ 18.5" ☐ 19.0" ☐ 19.5" ★ Not Available w/ 19.0", 19.5", 20.0", 20.5", 21.0" Front Seat Height *Measurement must be 1.5" shorter than your Front Seat Height	\$ STD														
HIGH MOUNT ☐ 10.0" ☐ 10.5" ☐ 11.0" ☐ 11.5" ☐ 12.0" ☐ 12.5" ☐ 13.0" ☐ 13.5" ☐ 14.0" ☐ 14.5" ☐ 15.0" ☐ 15.5" ☐ 16.0" With Risers or Clamps		\$ 65															
12. INSIDE WHEEL CLEARANCE		☐ 3/4 ☐ 1" ☐ 1 1/4" ☐ 1 1/2" ☐ 1 3/4" ☐ 2" *3/4" Inside wheel clearance not available with 2 degree camber. (with the exception of Y-360 wheels) *3/4" Inside wheel clearance not available with Armrest *3/4" Inside wheel clearance not available with D's Locks	\$ STD														
NOTE Razor blade is a custom made wheelchair. If you would like to customize your Razor Blade beyond what is available on this order form please submit your request here. Colours will determine as to whether this request can be accomplished while maintaining the safety and integrity of the wheelchair. If it can be accomplished the quote will reflect any and all custom charges, thank you.																	
Sub Total Section B 1/2			\$														
Sub Total Section B 2/2			\$														
Total Section B			\$														




C WHEELS, TIRES & CASTERS, FORKS 1/2

WHEELS														
Standard Radial Spoke Wheel ☐ 24" (540) Only		\$ STD												
Heat-Treated Wheel ☐ 24" (540) ☐ 25" (559) ☐ 26" (590)		\$ 250												
Cross Knot Laced Spokes ☐ 24" (540) ☐ 25" (559) ☐ 26" (590)		\$ 250												
Twister Super Strength Wheel ☐ 24" (540) ☐ 25" (559) ☐ 26" (590)		\$ 300												
SUN Fusion 16 Wheel ☐ 24" (540) ☐ 25" (559)		\$ 410												
16 Black Spokes Hub Color ☐ Black ☐ Silver														
Y-360 Billet Wheels ☐ 24" (540) ☐ 25" (559) *Max. Capacity 220Lbs *Camber 0- 2° max *Not for Sports or Rugged use Style ☐ Saw Blade ☐ Hearts(25"only) ☐ T-Style(25"only)		\$ 999												
Y-360 Billet Wheels MOD Style ☐ 25" (559) Only *Max. Capacity 220Lbs *Camber 0- 2° max *Not for Sports or Rugged use Style ☐ Tribal Sun ☐ Girls ☐ Eagle ☐ Holy Cross ☐ Hono Specify Powder Coat Color for Insert Only () *Choose from Frame Colors		\$ 1,499												
Extra Options for Y-360 ☐ Spinners Clover (\$650 up charge) ☐ Powder-Coated Spokes Only () *Choose from Frame Colors (\$150 up charge)														
X-Core (3 Spokes comp) ☐ 25" (559) Only *Camber 0- 2° max *Not available with D's Lock Wheel Color ☐ Black		\$ 743												
☐ Yellow ☐ Blue ☐ Red ☐ White		\$ 943												
Spinergy Spox ☐ 24" (540) ☐ 25" (559) ☐ 26" (590) 18 Spokes Spoke Color ☐ Black ☐ Red ☐ Yellow ☐ Blue ☐ White With Silver Hub With Black Hub With Black Hub With Black Hub With Black Hub		\$ 695												
☐ Custom Color Spoke Pattern <table><tr><td></td><td>Left Wheel</td><td>Right Wheel</td></tr><tr><td>Spokes</td><td></td><td></td></tr><tr><td>Color</td><td></td><td></td></tr><tr><td>Hub Color</td><td>☐ Black/☐ Silver</td><td>☐ Black/☐ Silver</td></tr></table>			Left Wheel	Right Wheel	Spokes			Color			Hub Color	☐ Black/ ☐ Silver	☐ Black/ ☐ Silver	\$ 795
	Left Wheel	Right Wheel												
Spokes														
Color														
Hub Color	☐ Black/ ☐ Silver	☐ Black/ ☐ Silver												
Spinergy LX ☐ 24" (540) ☐ 25" (559) ☐ 26" (590) 12 Spokes Spoke Color ☐ Black ☐ Red ☐ Yellow ☐ Blue With Silver Hub With Black Hub With Black Hub With Black Hub		\$ 895												
☐ White ☐ Neon Green ☐ Pink ☐ Orange With Black Hub With Black Hub With Black Hub With Black Hub														
☐ Custom Color Spoke Pattern <table><tr><td></td><td>Left Wheel</td><td>Right Wheel</td></tr><tr><td>Spokes</td><td></td><td></td></tr><tr><td>Color</td><td></td><td></td></tr><tr><td>Hub Color</td><td>☐ Black/☐ Silver</td><td>☐ Black/☐ Silver</td></tr></table>			Left Wheel	Right Wheel	Spokes			Color			Hub Color	☐ Black/ ☐ Silver	☐ Black/ ☐ Silver	\$ 995
	Left Wheel	Right Wheel												
Spokes														
Color														
Hub Color	☐ Black/ ☐ Silver	☐ Black/ ☐ Silver												
Topolino Cabon Core WX2.5 Wheels Black Spokes Only ☐ 24" (540) ☐ 25" (559) *Max. Capacity 220Lbs		\$ 1,095												
*Not available with D's Lock Nipples Color ☐ Red ☐ Yellow														
Pro Wheels ☐ 24" (540) ☐ 25" (559)		\$ 825												
w/ Aluminum Color Anodized Handrim Spoke & Handrim Color ☐ Anodized Black ☐ Anodized Red ☐ Anodized Blue *Not available with D's Lock														
HANDRIMS														
Type ☐ Standard For 24", 25" or 26" Wheel Anodized Black (*Not available with Pro Wheels)		\$ STD												
☐ Vinyl Coated For 24", 25", 26" Wheel (*Not available with Pro Wheels)		\$ 142												
☐ Natural Fit For 24", 25" Wheel ☐ Standard Grip (Not available with short tabs. Not available with Y-360 Billet Wheels, X-Core, Topolino Wheels, Pro Wheels, Mag Wheels)		\$ 295												
☐ Super Grip (Not available with short tabs. Not available with Y-360 Billet Wheels, X-Core, Topolino Wheels, Pro Wheels, Mag Wheels)		\$ 340												
☐ Without Thumb Grip (Not available with Y-360 Billet Wheels, Not available with Pro Wheels, Mag Wheels)		\$ 270												
☐ Natural Fit (LT) For 24", 25" Wheel ☐ Standard Grip (Not available with short tabs. Not available with Y-360 Billet Wheels, X-Core, Topolino Wheels, Pro Wheels, Mag Wheels)		\$ 325												
☐ Super Grip (Not available with short tabs. Not available with Y-360 Billet Wheels, X-Core, Topolino Wheels, Pro Wheels, Mag Wheels)		\$ 370												
☐ Without Thumb Grip (Not available with Y-360 Billet Wheels, Not available with Pro Wheels, Mag Wheels)		\$ 300												
☐ The Surge Handrims * Not available with Y-360 Billet Wheels * w / X-Core Wheels required Installation fee (\$30) flat fee		\$ 315												
☐ The Surge Handrims (LT) * Not available with Y-360 Billet Wheels * w / X-Core Wheels required Installation fee (\$30) flat fee		\$ 340												
☐ Spinergy Flex Rim Available only on 24" and 25" Spinergy LX.		\$ 495												
☐ Projections Handrims *Not available with Y-360 Billet Wheels, X-Core, Topolino Wheels, Pro Wheels, Mag Wheels) *Call For Availability		\$ CFP												
TAB LENGTH ☐ Long Tab (STD inside clearance 3/4") * Not available with X-Core Wheels		\$ STD												
☐ Short Tab (Inside clearance 1/2") *Not available with Ultra Grip and Natural Fit Handrim with thumb grip		\$ STD												
☐ Short Tab (Set closest to wheel cut off excess) * Not available with Y-360 wheel, Pro Wheels, X-Core, Pro Wheels		\$ 15												

C WHEELS, TIRES & CASTERS, FORKS 2/2

TIRES					
☐ High Pressure Tire	☐ 24"x1", ☐ 25"x1" or ☐ 26"x1"	Gray Tire with black wall Only		E2211/E2212	\$ STD
☐ High Pressure Tire	☐ 24"x1", ☐ 25"x1" or ☐ 26"x1"	☐ Blue ☐ Red		E2211/E2212	\$ 65
☐ Low Pressure Tire	☐ 24"x1-3/8" or ☐ 26"x1-3/8"	Gray Only		E2211/E2212	\$ STD
☐ E2211Schwalbe Marathon Plus High Pressure Tire	☐ 24"x1", ☐ 25"x1" or ☐ 26"x1"	Black Only		E2211/E2212	\$ 99
☐ Solid Tire	☐ 24"x1", ☐ 25"x1" or ☐ 26"x1"	☐ Dark Gray ☐ Black(\$30up charge)		E2218	\$ 135
☐ E2211 Moderate V-Track Knobby Tire	24"x1-3/8" (37-540)	Gray Only		E2211/E2212	\$ 117
☐ E2211 Aggressive Knobby Tire	24"x 2" (50-540)	Gray Only		E2211/E2212	\$ 117
☐ E2211 Aggressive Knobby Tire	26"x 1.95" (50-559)	Black Only		E2211/E2212	\$ 117
☐ E2213 Flat-Free Inserts For 24"x1-3/8" Tires Only				E2211/E2213	\$ 117
CASTERS					
SIZE	TIRE WIDTH	HUB COLOR	TIRE COLOR	DESCRIPTION	
☐ 3"	1"	Black Hub	Black Tire Only	Polyeurethane	\$ STD
☐ 3"	1"	Plastic Black Hub	Clear Tire Only	Lighted, Polyeurethane	\$ 65
☐ 3"	1"	Red Hub	Red Tire Only	Polyeurethane	\$ STD
☐ 3"	1"	Black Hub	Clear Tire Only	Polyeurethane	\$ STD
☐ 3"	1"	Aluminum Silver Hub	Black Tire Only	Made By Flog Legs	\$ 165
☐ 3"	1"	Aluminum Black Hub	Black Tire Only	Made By Flog Legs	\$ 165
☐ 3"	1.5"	Aluminum Silver Hub	Black Tire Only	Soft Roll Made By Flog Legs	\$ 165
☐ 3"	1.5"	Aluminum Black Hub	Black Tire Only	Soft Roll Made By Flog Legs	\$ 165
☐ 4"	1"	Plastic Black Hub	Black Tire Only		\$ STD
☐ 4"	1"	Plastic Black Hub	Clear Only	Lighted, Polyeurethane	\$ 65
☐ 4"	1"	Plastic Black Hub	Black Tire Only	Made By Flog Legs	\$ 125
☐ 4"	1"	Aluminum Silver Hub	Black Tire Only	Made By Flog Legs	\$ 165
☐ 4"	1"	Aluminum Black Hub	Black Tire Only	Made By Flog Legs	\$ 165
☐ 4"	1.5"	Composite Black Hub	Tire ☐ Black ☐ Blue ☐ Red	Soft Roll Made By Flog Legs	\$ 135
☐ 4"	1.5"	Aluminum Silver Hub	Tire ☐ Black ☐ Blue ☐ Red	Soft Roll Made By Flog Legs	\$ 165
☐ 4"	1.5"	Aluminum Black Hub	Black Tire Only	Soft Roll Made By Flog Legs	\$ 165
☐ 5"	1"	Composite Black Hub	Black Tire Only		\$ STD
☐ 5"	1"	Plastic Black Hub	Clear Tire Only	Lighted, Polyeurethane	\$ 65
☐ 5"	1"	Composite Black Hub	Black Tire Only	Made By Flog Legs	\$ 125
☐ 5"	1"	Aluminum Silver Hub	Black Tire Only	Made By Flog Legs	\$ 165
☐ 5"	1"	Aluminum Black Hub	Black Tire Only	Made By Flog Legs	\$ 165
☐ 5"	1.5"	Composite Black Hub	Tire ☐ Black ☐ Blue ☐ Red	Soft Roll Made By Flog Legs	\$ 135
☐ 5"	1.5"	Aluminum Silver Hub	Tire ☐ Black ☐ Blue ☐ Red	Soft Roll Made By Flog Legs	\$ 165
☐ 5"	1.5"	Aluminum Black Hub	Black Tire Only	Soft Roll Made By Flog Legs	\$ 165
FORKS					
☐ STD Aluminum Forks	Black Only				\$ STD
Extra Options for Aluminum Forks	☐ With Suspension Springs		E1015	\$ 55	
	☐ With Quick Release Axles		K0073	\$ 107	
☐ Frog Leg Suspension Forks Std 1" Wide (For 3", 4" x 1" width Caster and 5" STD Caster, 5" lighted Caster)	Color ☐ Anodized Black ☐ Anodized Silver		E1015	\$ 575	(Not available with Quick Release Axles)
☐ Frog Leg Suspension Forks 1.5" Wide (For 3", 4" x 1.5" width Caster and 5" x 1" & 5" x 1.5" Casters)	Color ☐ Anodized Black ☐ Anodized Silver		E1015	\$ 575	(Not available with Quick Release Axles)
Sub Total Section C 1/2				\$	
Sub Total Section C 2/2				\$	
Total Section C				\$	

D COMPONENTS

WHEEL LOCKS								
☐ Push to Lock / Anodized Black Only		\$ STD						
☐ NEW Compact Push to Lock / Anodized Black Only ☐ w/ 2" Extension Handle (\$65up)		\$ 110						
☐ NEW Compact Scissor Lock / Anodized Black Only		\$ 110						
☐ D's Locks (*Not available with X-Core, Topolino WX2.5, Pro Wheels, Magic wheels) (*Does not include installation fee)		\$ CFP						
PUSH HANDLES *If you are ordering push handles, You must select the backrest locking mechanism.								
☐ Integrated w/Fixed Backrest Height / Anodized Black Only		\$ STD						
☐ Bolt-on / Anodized Black Only Specify Height from floor -()		\$ 142						
☐ Bolt-on w/Quick Release / Anodized Black Only Specify Height from floor -()		\$ 190						
☐ NEW Height Adj. / Anodized Black Only Specify Maximum height from floor -()		\$ 295						
NEW ARM RESTS Anodized Black Only,								
☐ Fixed Tubular Armrest	Adj. Angle and Height Height from seat sling Set at : ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" E0973	\$ 250						
☐ Swing Away Tubular Armrest	Adj. Angle and Height Height from seat sling Set at : ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" E0973	\$ 250						
☐ Flip-Up Armrest	☐ w/Fixed Tubular Armrest Height from seat sling Set at : ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" E0973	\$ 250						
Adj. Angle	☐ w/ Swing Away Tubular Armrest Height from seat sling Set at : ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" E0973	\$ 250						
☐ Custom Length & Height (\$150 up charge) *Specify Height & Length → 								
ANTI-TIPPERS								
☐ Anti-Tippers (pair) w/ Quick Release / Anodized Black Only		E0971 \$ 112						
SIDE GUARDS								
☐ Rigid Side Guards - Black ABS	☐ Non-Removable ☐ Removable ☐ Removable w/Hole (\$10 up charge) Size: ☐ STD ☐ Medium ☐ Cut to Wheel (\$98 up charge) Measurement Color: Black Only	\$ 142						
☐ Rigid Side Guards - Aluminum  <table border="1"> <tr> <th>Height</th> <th>Size</th> </tr> <tr> <td>7.5"</td> <td>STD</td> </tr> <tr> <td>8.25"</td> <td>Medium</td> </tr> </table>	Height	Size	7.5"	STD	8.25"	Medium	☐ Non-Removable ☐ Removable Size: ☐ Medium Only ☐ Cut to Wheel (\$98 up charge) Measurement Color ☐ Black (STD) ☐ Pol. Aluminum (STD) ☐ Powder Coated () (\$50 up charge) ☐ Gold Plate (+CFP)	\$ 197
Height	Size							
7.5"	STD							
8.25"	Medium							
								
☐ PRO-TEX Fender Guards / (Black ABS -Only) ☐ Cut to Fit Specs ☐ Do Not Install		\$ 175						
☐ Fabric Clothing Guards Black Nylon w/Velcro		\$ 78						
SEAT CUSHION *Please see section E for other cushion choices.-								
☐ Foam Cushion w/Nylon Black Cover Thickness ☐ 1" ☐ 2" ☐ 3" ☐ 4"		\$ 78						
☐ Custom Cushion Cover (Std. Foam Cushion Only) -Select Fabric () (\$90 up charge)								
MISCELLANEOUS								
☐ Frame Protectors / Black Leather/Zipper Type		\$ 78						
☐ K0065 Spoke Protectors / Black Only ☐ 24" ☐ 25" ☐ 26"		K0065 \$ 107						
☐ E0978 Seat Belt w/Auto Buckle (BP) 1.5" Black Belt Only		E0978 \$ 69						
☐ Tie Down Brackets 4pcs.		\$ 210						
SPECIALTY PARTS ☐ Wheelchair Tray Size: ☐ W14"x D12" ☐ W16"x D14" ☐ W18"x D14" Color: ☐ Raw Alum. ☐ Black Alum. ☐ Polish Alum. (\$250 up charge)								
\$ 250								
UPHOLSTERY CHOICE & COLOR *Positioning Back Available. Please see Section E								
☐ Standard Adjustable Ht Backrest Black Only		\$ STD						
☐ Adjustable Tension Backrest Black Only (*Not available with 8" Backrest Height)		\$ 135						
Custom Options	☐ Fabrics for Fun ☐ 001 ☐ 002 ☐ 003 ☐ 004 ☐ 006 ☐ 010 ☐ 011 . . . (\$95 up charge) For Std. or Adj. Style ☐ 014 ☐ 015 ☐ 016 ☐ 017 ☐ 018 ☐ 020 ☐ 021 ☐ 022 ☐ 023 ☐ 024 ☐ 025 ☐ 026 ☐ 027 ☐ 028 ☐ Backrest Leather Trim (outside panels) (\$77 up charge)							
☐ Calf Strap ☐ Black ☐ Matching Fabric (\$55 up charge) -Select Fabric()		K0038 \$ 24						
☐ Matching Fabric on Padded Side Guards -Select Fabric ()		\$ 210						
Total Section D \$								

CHAIR OPTIONS 1/2

SEAT CUSHION									
SUPRACOR® WHEELCHAIR CUSHIONS									
Stimulite® Contoured Cushion		☐ W14"xD16" (CD1416)						E2603	\$ 525
*Black Upholstery only		☐ W16"xD16" (CD1616)						E2603	\$ 525
		☐ W16"xD18" (CD1618)						E2603	\$ 525
		☐ W18"xD16" (CD1816)						E2603	\$ 525
		☐ W18"xD18" (CD1818)						E2603	\$ 525
		☐ W18"xD20" (CD1820)						E2603	\$ 540
		☐ W20"xD16" (CD2016)						E2603	\$ 540
		☐ W20"xD18" (CD2018)						E2603	\$ 540
		☐ W20"xD20" (CD2020)						E2603	\$ 555
		☐ Extra Sensitive (\$30 up charge)							
		☐ Sling Style Bottom (\$45 up charge)							
		☐ Custom Size W(") x D(") (\$- CFP)							
Stimulite® Slimline Cushion		☐ W14"xD16" (SM1416)						E2603	\$ 495
*Black Upholstery only		☐ W16"xD16" (SM1616)						E2603	\$ 495
		☐ W16"xD18" (SM1618)						E2603	\$ 495
		☐ W18"xD18" (SM1818)						E2603	\$ 495
		☐ W18"xD20" (SM1820)						E2603	\$ 495
		☐ W20"xD16" (SM2016)						E2603	\$ 510
		☐ W20"xD18" (SM2018)						E2603	\$ 510
		☐ W20"xD20" (SM2020)						E2603	\$ 510
		☐ Extra Sensitive (\$30 up charge)							
		☐ Sling Style Bottom (\$45 up charge)							
		☐ Custom Size W(") x D(") (\$- CFP)							
Stimulite® Corbee Pediatric		☐ W10"xD10" (PD1010)						E2603	\$ 450
*Black Upholstery only		☐ W12"xD12" (PD1212)						E2603	\$ 450
		☐ W14"xD14" (PD1414)						E2603	\$ 450
		☐ Extra Sensitive (\$30 up charge)							
		☐ Sling Style Bottom (\$45 up charge)							
		☐ Custom Size W(") x D(") (\$- CFP)							
ROHO WHEELCHAIR CUSHIONS									
4" HIGH PROFILE® QUADTRO SELECT® CUSHION		*Black Upholstery						E2622-2623	\$ 525
Choose Size	☐ W12"xD12"	☐ W12"xD13"	☐ W13"xD13"	☐ W13"xD15"	☐ W13"xD16"	☐ W15"xD13"	☐ W15"xD15"		
	☐ W15"xD16"	☐ W15"xD18"	☐ W16"xD15"	☐ W16"xD16"	☐ W16"xD18"	☐ W16"xD20"	☐ W18"xD15"		
	☐ W18"xD16"	☐ W18"xD18"	☐ W18"xD20"	☐ W18"xD22"	☐ W20"xD16"	☐ W20"xD18"	☐ W20"xD20"		
	☐ W20"xD22"	☐ W22"xD18"	☐ W22"xD20"	☐ W22"xD22"	☐ Custom Size W(") x D(") (\$+85)				
3" MID PROFILE® QUADTRO SELECT® CUSHION		*Black Upholstery only						E2622-2623	\$ 525
Choose Size	☐ W15"xD15"	☐ W15"xD16"	☐ W16"xD16"	☐ W16"xD18"	☐ W18"xD16"	☐ W18"xD18"			
	☐ Custom Size W(") x D(") (\$+85)								
2" LOW PROFILE® QUADTRO SELECT® CUSHION		*Black Upholstery only						E2622-2623	\$ 525
Choose Size	☐ W12"xD12"	☐ W12"xD13"	☐ W13"xD13"	☐ W13"xD15"	☐ W13"xD16"	☐ W15"xD13"	☐ W15"xD15"		
	☐ W15"xD16"	☐ W15"xD18"	☐ W16"xD15"	☐ W16"xD16"	☐ W16"xD18"	☐ W16"xD20"	☐ W18"xD15"		
	☐ W18"xD16"	☐ W18"xD18"	☐ W18"xD20"	☐ W20"xD16"	☐ W20"xD18"	☐ W20"xD20"			
	☐ Custom Size W(") x D(") (\$+85)								
2" NEXUS SPIRIT® CUSHION		*Black Upholstery only						E2622	\$ 525
Choose Size	☐ W14"xD14"	☐ W14"xD16"	☐ W15"xD15"	☐ W16"xD16"	☐ W16"xD18"	☐ W17"xD17"	☐ W18"xD16"		
	☐ W18"xD18"	☐ W20"xD18"	* Custom Size Not Available						
STAR CUSHIONS									
Galaxy Cushion *Black Upholstery only						Choose Zoning Options		E2624	\$ 498
☐ W14"xD14" (GXY 1414-2)		☐ W14"xD16" (GXY 1416-2)		⇒	☐ Standard - Dual Valve				
☐ W16"xD14" (GXY 1614-2)		☐ W16"xD16" (GXY 1616-2)			☐ Dual Valve - Front to Back				
☐ W16"xD18" (GXY 1618-2)		☐ W18"xD16" (GXY 1816-2)			☐ Four Valve				
☐ W18"xD18" (GXY 1818-2)		☐ W18"xD20" (GXY 1820-2)			☐ Simple - Single Valve				
☐ W20"xD16" (GXY 2016-2)		☐ W20"xD18" (GXY 2018-2)			☐ Four Valve - Enhanced				
☐ W20"xD20" (GXY 2020-2)					☐ Four Valve - Front to Back				
THE CURVE CUSHION		☐ W10"xD10"	☐ W10"xD11"	☐ W10"xD11"	☐ W10"xD12"	☐ W11"xD12"	☐ W12"xD12"	\$ 133	
*Black Upholstery only		☐ W12"xD13"	☐ W12"xD14"	☐ W13"xD11"	☐ W13"xD12"	☐ W13"xD13"	☐ W13"xD14"		
		☐ W14"xD14"	☐ W14"xD15"	☐ W14"xD16"	☐ W15"xD14"	☐ W15"xD15"	☐ W15"xD16"		
		☐ W16"xD14"	☐ W16"xD15"	☐ W16"xD16"	☐ W16"xD17"	☐ W16"xD18"	☐ W17"xD16"		
		☐ W17"xD17"	☐ W17"xD18"	☐ W18"xD16"	☐ W18"xD17"	☐ W18"xD18"			
NOTE									

E CHAIR OPTIONS 2/2

POSITIONING BACK			
*Note : Actual height of the backrest will increase 2" or more.			
ADI CARBON BACKREST (Support Contour-ACTIVE 2.5") *Not Available with Welded Side Guards Mounting Hardware: ☐ 2 Point Pro ☐ Quick Release			
Height	☐ LOW (10" Tall) *Backrest Height requires 12" or more *Available 13"-19" width ☐ Medium (13" Tall) *Backrest Height requires 15" or more *Available 13"-19" width ☐ Standard (16" Tall) *Backrest Height requires 18" or more *Available 13"-19" width ☐ Custom Size (" Tall) x (" width)	Width (") E2613	\$ 745
		Width (") E2613	\$ 815
		Width (") E2613	\$ 880
		E2613	\$ CFP
ADI ALUMINUM BACKREST (Support Contour-ACTIVE 2.5") *Not Available with Welded Side Guards Mounting Hardware: ☐ 2 Point Pro ☐ Quick Release			
Height	☐ LOW (10" Tall) *Backrest Height requires 12" or more *Available 10"-20" width ☐ Medium (13" Tall) *Backrest Height requires 15" or more *Available 12"-20" width ☐ Standard (16" Tall) *Backrest Height requires 18" or more *Available 14"-20" width ☐ Tall (18" Tall) *Backrest Height requires 20" or more *Available 16"-20" width ☐ X-Tall (20" Tall) *Backrest Height requires 22" or more *Available 16"-20" width ☐ Custom Size (" Tall) x (" width)	Width (") E2613	\$ 610
		Width (") E2613	\$ 625
		Width (") E2613	\$ 640
		Width (") E2613	\$ 655
		Width (") E2613	\$ 670
		E2613	\$ CFP
ADI ALUMINUM BACKREST (Support Contour-DEEP 4") *Not Available with Welded Side Guards *Not Available with Lateral Supports Mounting Hardware: ☐ 2 Point Pro ☐ Quick Release			
Height	☐ LOW (10" Tall) *Backrest Height requires 12" or more *Available 10"-20" width ☐ Medium (13" Tall) *Backrest Height requires 15" or more *Available 12"-20" width ☐ Standard (16" Tall) *Backrest Height requires 18" or more *Available 14"-20" width ☐ Tall (18" Tall) *Backrest Height requires 20" or more *Available 16"-20" width ☐ X-Tall (20" Tall) *Backrest Height requires 22" or more *Available 16"-20" width ☐ Custom Size (" Tall) x (" width)	Width (") E2615	\$ 710
		Width (") E2615	\$ 725
		Width (") E2615	\$ 740
		Width (") E2615	\$ 755
		Width (") E2615	\$ 770
		E2615	\$ CFP
ADI LATERAL SUPPORTS (Lateral Supports are compatible with only active 2.5" ADI positioning back.) ☐ Stationary Pad Pad Size ☐ Small (3" H x 3" D) ☐ Medium(3" H x 4" D) ☐ Large(4.5" H x 4.5" D) ☐ Large(4" H x 6" D) ☐ Pair \$ 320 ☐ Left side only or ☐ Right side only \$ 160 ☐ Swing-Out Pad Size ☐ Small (3" H x 3" D) ☐ Medium(3" H x 4" D) ☐ Large(4.5" H x 4.5" D) ☐ Large(4" H x 6" D) ☐ Pair \$ 520 ☐ Left side only or ☐ Right side only \$ 260			
OTHER			
☐ ADI Transfer Board ☐ 21" (ADV21) ☐ 24" (ADV24) ☐ 29"			\$ 62
☐ ADI Push and Transfer Gloves Leather (one size fits)			\$ 60
☐ ADI Carbon Solid Seat (Wide: 18" Max, Depth: 20") ※Colours Side Guards Not Available with ADI Carbon Solid Seat			\$ 625
With / ☐ ADI Folding Carbon Side Guard			\$ 400
With / ☐ ADI Fixed Carbon Side Guard			\$ 375
VARILITE® Seating Systems Available			CFP
Please visit for all options. www.varilite.com			
NOTE			
Sub Total Section E 1/2 \$			
Sub Total Section E 2/2 \$			
Total Section E \$			

F SPECIAL PARTS COLORS (Please allow an additional 1-2 weeks delivery)

PARTS COLORS

Powder-Coated	Gloss Black	Twilight Flake	Yellow	Magenta	Fire Engine Red	White
	Candy Blue	Candy Purple	Candy Green	Candy Red	Candy Raspberry	Cobalt Blue
	Matte Gray	Matte Purple	Matte Blue	Matte Green	Matte Black	Matte Red
	Neon Red	Neon Green	Orange Sparkle	Mellon Yellow	Gun Smoke Silver	Pigeon Blue
	Granny Sparkle Green	Sparkle Sky Blue	Vermilion	Reed Green	Blue Lilac	

Special Colors Polished Aluminum
24 karat Gold Plated

Please Choose the Color for each Part

*Parts not selected will be black color. (Suspension Spring: Silver Color)

*Minimum 3 parts selected per order.

*PC=Powder-Coated

Saddle Stops

- ☐ PC () \$40
- ☐ Polished..... \$115
- ☐ Gold Plated..... \$CFP

Backrest

- ☐ PC () \$125
- ☐ Polished..... \$225
- ☐ Gold Plated..... CFP

Backrest Brackets

- ☐ PC () \$90
- ☐ Polished..... \$115
- ☐ Gold Plated..... \$CFP

Foot Plate

- ☐ PC () \$45
- ☐ Polished..... \$90
- ☐ Gold Plated..... \$CFP

Footrest Rings

- ☐ PC () \$20
- ☐ Polished..... \$40
- ☐ Gold Plated..... \$CFP

Foot Plate Clamps

- ☐ PC () \$50
- ☐ Polished..... \$75
- ☐ Gold Plated..... \$CFP

Footrest Tubes

- ☐ PC () \$65
- ☐ Polished..... \$85
- ☐ Gold Plated..... \$CFP

Caps

- ☐ PC () \$20
- ☐ Polished..... N/A
- ☐ Gold Plated..... N/A

Front Springs

- ☐ PC () \$20
- ☐ Polished..... N/A
- ☐ Gold Plated..... N/A

Forks & Knurl Nuts

- ☐ PC () \$65
- ☐ Polished..... \$90
- ☐ Gold Plated..... \$CFP

Axle Plates

- ☐ PC () \$65
- ☐ Polished..... \$85
- ☐ Gold Plated..... \$CFP

Camber Plates

- ☐ PC () \$50
- ☐ Polished..... \$75
- ☐ Gold Plated..... \$CFP

Axle Bases

- ☐ PC () \$50
- ☐ Polished..... \$75
- ☐ Gold Plated..... \$CFP

Axle Clamp

- ☐ PC () \$50
- ☐ Polished..... \$75
- ☐ Gold Plated..... \$CFP


Axle Tube

- ☐ PC () \$50
- ☐ Polished..... \$75
- ☐ Gold Plated..... \$CFP


Total Section F \$

1. SEAT WIDTH


Measure from outside of seat tube to the outside of the opposite seat tube.

**2. SEAT DEPTH**


Measure from at the front bottom of backrest tube to the front edge of seat sling.

**3. FOOTREST WIDTH**


Measure from inside of front frame tube to inside of opposite front frame tube.

**4. FRAME ANGLE**


Measure angle frame front to floor.

**5. BACKREST HEIGHT**


Measure from top of back post to seat sling at front bottom of backrest tube.

**6. BACKREST ANGLE**


Measure angle from front of back post to floor.

**7. FRONT SEAT HEIGHT**

Measure from floor to top of seat sling at front edge.


**8. REAR SEAT HEIGHT**

Measure from floor to top of seat sling at front bottom of backrest tube.


**9. CENTER OF GRAVITY**

Measure from center of back tube at the folding pivot to the center of the rear axle.


(If welding backrest is selected, Measure from center of back tube at the seat tube to the center of the rear axle.)

**10. CAMBER ANGLE**

Measure angle from wheel to vertical at floor.

**11. FRONT OF SEAT SLING TO FOOTREST****ADULT CHAIR**

Measure from top of seat sling at front edge to top of foot plate at front edge.

**PEDIATRIC CHAIR**

Measure vertical from top of seat sling at front edge to top of foot plate.

**12. INSIDE WHEEL CLEARANCE**

Measure from outside of seat back tube to inside of rear tire.

(If side guards are selected, measure from outside of side guards to inside of rear tire.)

