

XO-202

Power-Standing/Drive Lightweight Wheelchair



- **USA Patent /CE Approved**
- **6061 T-6 Aircraft-grade Aluminum Frame**
- **Flip-back Padded Armrests w/SHARK II Controller & Swivel Functions For Easy Access + Shorter Overall Length**
- **High Powered Lloyd Control Box**
- **On/Off Master Switch**
- **Adjustable Backrest Height**
- **Adjustable Footrest Height**
- **Adjustable Knee Support, Safety Belt, Chest Support & Leg Strap**
- **2-3" Memory Foam Cushion Back and Seat**
- **Weight Only: 110 lbs w/o Battery(lbs w/ Battery)**
- **Breathable Nylon Upholstery Seat And Back w/ Pouch/Zipper Fanny Pack**
- **See Through Sleek Side Panels**
- **Flat Free Front & Rear Wheels**
- **Motor: DC24V. 320W**
- **Maximum Speed: 5 Mph (Moving slowly is recommended while standing up)**
- **14" High Strength Flat Free Rear Tires**
- **8" Flat Free Front Caster**
- **Battery: 12V/36AH x2pcs**
- **Range per Charge: ~25 Miles**(depending on weight of user and terrain)

**Options: Push Handles, Adjustable Headrest, and Vehicle Tie Downs Kit

Specifications:

Available In A Lighter Manual Propulsion Model (XO-101)!

XO-202N → **Seat Width (16")**, **Seat Depth (18/19/20")**, **Seat Height (20 ½")**, **Back Height (16-19")**, **Overall Height (35 - 56 ½")**, **Overall Width (24 ½")**, **Overall Length (40")**, **Armrest Height (6 ¼" – 10 ¾")**, **Turning Radius (25")**, **Weight Capacity (250 lbs)**

XO-202 → **Seat Width (18")**, **Seat Depth (18/19/20")**, **Seat Height (20 ½")**, **Back Height (16-19")**, **Overall Height (35 – 56 ½")**, **Overall Width (26 ½")**, **Overall Length (40")**, **Armrest Height (6 ¼" – 10 ¾")**, **Turning Radius (25")**, **Weight Capacity (250 lbs)**

*Karman Healthcare Inc. Recommends Patient
Seek Physicians Counsel Before Purchase



Karman Healthcare Inc.
19255 San Jose Ave.
City Of Industry, CA 91748-1418
Order: 1-800-805-2762
Tel: 626-581-2235
Fax: 626-581-2335
www.karmanhealthcare.com
Email: info001@karmanhealthcare.com