



S-Ergo 115-TP

Ultra Lightweight Ergonomic Manual Transporter

- S-Style Ergonomic Seat
- High strength, weighs only 22 lbs. (w/o footrests)
- Silver ¼” Removable Anti-bacterial cushion
- Folding backrest/folding seat for transporting in vehicle or travel
- Fixed armrests w/ wider concaved armpads
- Swing In & Away footrests allowing maximum safety while entering or exiting the chair.
- “Tube-in Center” foot-plate, assures better side leg support
- Extending push to lock manual brakes
- Pouch for carrying small items attached to upholstery behind backrest.
- 14” Rear polyurethane, high profile, flat free wheels
- 7”x1” Polyurethane front caster
- Available in 16”, 18” or 20” seat width
- Frame Color: Pearl Silver & Rose Red
- Upholstery: Black Breathable Mesh Bottom
- Weight Capacity of 250 lbs

Features for S-Ergo Series

- Ergonomic, ultra lightweight specialized frame
- Pressure relief, stabilization, and anti-slipage seating system
- Upholstery is equipped with AEGIS anti-biotic technology, odor/stain/deterioration control.
- Wider/concaved armpad, providing comfort & support

	seat width	seat depth	armrest height	seat height	back height	overall height	overall open width	folded width	overall length	weight w/o riggings	weight capacity
Model No.	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(lb.)	(lb.)
S-115TP16	16	17	8	19	16	36	22	12	36	22	250
S-115TP18	18	17	8	19	16	36	23	12	36	22	250
S-115TP20	18	17	8	19	16	36	23	12	36	23	250

*Available Options ~ Anti-Tipper, Seat belt, Extra Upholstery, Memory Foam, UV Resistant Canopy, Rigidifying Push-Bar, Reflecting Spoke Guard, Foldable/Adjustable Headrest, Frog Leg Suspension, Magnesium Wheel, Vulcan Grip Wheel

* S-115TP-WB – Wire-Brake Push Handle model available in 16” or 18” seat width only.



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue. City of Industry. CA 91748-1418
 Order. 1-800-805-2762 Tel. 626-581-2235 Fax. 626-581-2335