





N472-MU


①

②





③

④

⑤

⑥

⑦



①	Back uph ass'y	4370CM40(26" BLACK NYLON SEAT UPHOLSTERY)
	4370CM43(26" BLACK NYLON BACK UPHOLSTERY)	4370CM41(28" BLACK NYLON SEAT UPHOLSTERY)
	4370CM44(28" BLACK NYLON BACK UPHOLSTERY)	4370CM42(30" BLACK NYLON SEAT UPHOLSTERY)
	4370CM45(30" BLACK NYLON BACK UPHOLSTERY)	
②	Rear wheel Ass'y.	⑦ Swingaway footrest ass'y
	3250AM44(24"*1-3/8"*ST REAR WHEEL ASS'Y 700LBS).	3330AM74(CP SWINGANG FOOTREST ASS'Y LH/RH)
③	Brake ass'y	3330AM75(CP SWINGANG FOOTREST ASS'Y LH)
	33200M67 BRAKE ASS'Y LH	3330AM76(CP SWINGANG FOOTREST ASS'Y RH)
	33200M68 BRAKE ASS'Y RH	33304M72(BLACK SWINGANG FOOTREST ASS'Y LH/RH)
④	Front wheel ass'y	33304M79(BLACK SWINGANG FOOTREST ASS'Y LH)
	33001M09(FRONT WHEEL ASS'Y)	33304M80(BLACK SWINGANG FOOTREST ASS'Y RH)
⑤	Detachable desk arm ass'y LH/RH	
	3100BM32(DETACHABLE DESK ARM ASS'Y LH/RH)	
	3100BM33(DETACHABLE DESK ARM ASS'Y LH)	
	3100BM34(DETACHABLE DESK ARM ASS'Y RH)	
⑥	Seat uph ass'y	