

P310--ARMU/P310--SBMU/P310--SDMU



①	Seat ass'y
	32200060 (Seat ass'y)
	32222011 (CAPTAIN SEAT ASS'Y,18W*16D*20H(M12B))
	43701072(PAD, ARMREST(LONG),BLACK)
	43701073(PAD,ARMREST(SHORT),BLACK)
②	Shroud ass'y
	30510030 (SHROUD SET, CANDY APPLERED, MP-3C)
	30515018 (SHROUD SET, MP-3C,COBALT BLUE)
	30552011 (SHROUD ASS'Y., SILVERBLACK, MP-3C)
③	Rear wheel ass'y.
	32600431(REAR WHEEL ASS'Y)
④	Motor and gearbox complete, LH and RH.
	62010074 (MOTOR & GEARBOX ASS'Y.,LH,P310)
	62020074 (MOTOR & GEARBOX ASS'Y.,RH,P310)

	62110065 (MOTOR,M1CG/H,L.)
	62120065 (MOTOR SET,M1CG/H,R.)
	62410027 (GEARBOX ONLY,LH,MP3C)
	62420027 (GEARBOX ONLY,RH, MP3C)
⑤	Controller
	33703023 (DYNAMIC SHARK CONTROLLER)
	64810005 (JOYSTICK ONLY, DY SHARK,DK-REMA03)
	64210005(CONTROLLER ONLY,DY SHARKDK-PMA05)
⑥	Caster Wheel with Fork Ass'y.
	33004002(200*50 PU CASTER WHEELASS'Y., LH/RH)
⑦	Footplate ass'y.
	33303037(ALUM. FOOTPLATE&BRACKETASS'Y.)
⑧	Charger unit
	65142005 (CHARGER, HP8204-5A-OFF)
	63600008 (CHARGER POWERCORD W/PLUG,US SPEC)