


P314--ARMU/P314--SBMU/P314--SPMU



①	Seat ass'y		33100609 (CASTOR ASS'Y.,200x50 PU WHEEL, RH)
	32200141 (Seat ass'y)	⑤	Controller
	32228001 (Seat Ass'y 20"X18"X20"(Mh82 Black Gray)		33503058 (Joystick Bracket Ass'y.,Shark, Dk-Rema01, ED)
	43701072(PAD, ARMREST(LONG),BLACK)		64810012 (JOYSTICK, DY SHARK,DK-REMA01)
	43701073(PAD,ARMREST(SHORT),BLACK)		64210005 (Controller Only,Dy Sharkdk-Pma05)
②	Drive Wheel Ass'y.	⑥	Shroud ass'y
	32600432 (DRIVER WHEEL ASS'Y.,260*85 PU)		30500067 (CENTER COVER ASS'Y BLACK, COLOR)
③	Motor and gearbox complete, LH and RH.		30500066 (REAR SHROUD, BLACK,P313/P314)
	62010116 (MOTOR & GEARBOX, LHM11LA-25)		30510347 (MAIN SHROUD, RED, P313)
	62020116 (MOTOR & GEARBOX (RH)M11LA-25)		30515180 (MAIN SHROUD, BLUE, P313)
	62110071 (MOTOR,M1LA/B,L.)		30521017 (MAIN SHROUD, CHAMPAGNE, P313)
	62120071 (MOTOR,M1LA/B,R.)	⑦	Footplate ass'y.
	62410059 (GEARBOX ASS'Y, M11-25,L.)		33300341 (FOOTPLATE ASS'Y.W/MERITS LOGO)
	62420059 (GEARBOX ASS'Y, M11-25,R.)	⑧	Charger unit
④	Caster Wheel with Fork Ass'y.		65142005 (CHARGER, HP8204-5A-OFF)
	33100607 (CASTOR ASS'Y.,200x50 PU WHEEL, LH)		63600008(CHARGER POWERCORD W/PLUG,US SPEC)