

# P32011ARMU/P32011SBMU



①	Seat ass'y 32200150 (Seat ass'y) 32222079 (SEAT ASS'Y.,MERITSFISH-ON DX, SB956) 43700401 (PAD, ARMREST,BLACK)	⑤	Caster Wheel with Fork Ass'y. 33100152 (CASTOR FORK ASS'Y., LH7*1-3/4" PU) 33100153 (CASTOR FORK ASS'Y., RH7*1-3/4" PU)
②	Battery Pack Ass'y 33700068(BATTERY PACK ASS'Y.,22AH, BLACK, P320)	⑥	Controller 33603032 (JOYSTICK BRAKET ASS'Y.,SHARK, ED) 64810006 (JOYSTICK, DY SHARK,DK-REMA02) 64210005(CONTROLLER ONLY,DY SHARKDK-PMA05)
③	Drive wheel ass'y 32600411 (DRIVER WHEEL, 2.80/2.50-4 PU,BLACK)	⑦	Shroud ass'y 30510313(SHROUD ASS'Y., CANDYAPPLE, P320) 30515159(SHROUD ASS'Y., BLUE,P320)
④	Motor and gearbox complete, LH and RH. 62011006 (MOTOR & GEARBOX ASS'Y.,M72AC-19, LH) 62021006 (MOTOR & GEARBOX ASS'Y.,M72AC-19, RH) 62111002 (MOTOR ASS'Y., M7AC/D, LH) 62121002 (MOTOR ASS'Y., M7AC/D, RH) 62411004 (GEARBOX ASS'Y., M73-19,LH ) 62421004 (GEARBOX ASS'Y., M73-19,RH)	⑧	Footplate ass'y. 42500157(FOOTREST PLATE)
		⑨	Charger unit 65142007(CHARGER, HP8204-3A-OFF) 63600008(CHARGER POWERCORD W/PLUG,US SPEC)