

P327--ARMU/P327--SBMU


①	Seat ass'y 32200175 (Seat ass'y) 32222085 (Seat Ass'y20.5W18D20H MC7/15/22/25) 43701072(PAD, ARMREST(LONG),BLACK) 43701073(PAD,ARMREST(SHORT),BLACK)		
②	Drive Wheel Ass'y. 32600423 (DRIVER WHEEL ASS'Y.,4.00-5 C-154 PT)		
③	Anti-tipper arm ass'y., Rear 30200247 (Anti-tipper arm ass'y., Rear, LH) 30200248 (Anti-tipper arm ass'y., Rear, RH)		30200246 (Anti-tipper arm ass'y., Front, LH)
④	Motor and gearbox complete, LH and RH. 62012007 (Motor & Gearbox ass'y., M9L01DA-5, LH) 62022007 (Motor & Gearbox ass'y., M9L01DA-5, RH) 62112007 (Motor ass'y., M9DA-5, LH) 62122007 (Motor ass'y., M9DA-5, RH) 62412001 (Gearbox ass'y., M9L01,L.) 62422001 (Gearbox ass'y., M9L01,R.)		⑥ Controller 33503064 (JOYSTICK ASS'Y.,DY Shark,SWING-OUT) 64810013 (JOYSTICK ONLY, DY Shark DK-REMD01) 64210010 (CONTROLLER ONLY, DY DK-PMB40).
⑤	Anti-tipper arm, Front 30200245 (Anti-tipper arm ass'y., Front, LH)		⑦ Shroud ass'y 44600508 (SHROUD, BLACK, P327, LH) 44600509 (SHROUD, BLACK, P327, RH) 43600608 (REAR SHROUD, BLACK, P327) 30510899 (SHROUD, FENDER, RED, P327, LH) 30510900 (SHROUD, FENDER, RED, P327, RH) 30515990 (SHROUD, FENDER, BLUE, P327, LH) 30515989 (SHROUD, FENDER, BLUE, P327, RH) 30500088 (MAIN SHROUD, BLACK, P327)
			⑧ Footplate ass'y. 33300350 (FOOTPLATE ASS'Y.W/MERITS LOGO)
			⑨ Charger unit 65142005 (CHARGER, HP8204-5A-OFF) 63600008(CHARGER POWERCORD W/PLUG,US SPEC)