

P7102-ARMU/P7102-SBMU



①	Seat ass'y	62412001(GEARBOX ASS'Y, M9L01,LH)
	32200179 (CAPTAIN SEAT ASSY)	62422001(GEARBOX ASS'Y, M9L01,RH)
	32222089 (CAPTAIN SEAT 22W*18D*22H)	⑤
	43701072(PAD, ARMREST(LONG),BLACK)	33703053(JOYSTICK ASS'Y,VR2,SWING-OUT)
	43701073(PAD,ARMREST(SHORT),BLACK)	64820002 (JOYSTICK ONLY, P&G VR2,JSM D50677)
②	Drive Wheel Ass'y.	64220003 (CONTROLLER ONLY, PG VR2-60, D50681)
	32600423(12" Rear Wheel Ass'y.,4.00-5 Alu. Rim, Black,Pt)	⑥
③	Caster Wheel with Fork Ass'y.	30510834 (SHROUD ASS'Y., CANDYAPPLE, P710)
	33100627(Castor Wheel Ass'y., 2.80/2.50-4 Puwheel)	30515191 (SHROUD ASS'Y., BLUE,P710)
④	Motor and gearbox complete, LH and RH.	⑦
	62012002 (MOTOR&GEARBOX ASS'Y,M9L01AB-1, LH)	33303059 (FOOTPLATE ASS'Y, P720)
	62022002 (MOTOR&GEARBOX ASS'Y,M9L01AB-1, RH)	⑧
	62112002 (MOTOR ASS'Y, M9AB-1,LH)	65102001(Charger, Off Board,Hp8204-6a-Off,115/230VAC)
	62122002 (MOTOR ASS'Y, M9AB-1,RH)	63600008(CHARGER POWERCORD W/PLUG,US SPEC)