

SUITCASE® Signature Series®



The SUITCASE® Signature Series® portable wheelchair ramp, with its single-fold design, sets up quickly and can be carried like a suitcase when not in use! It features a no-pinch, non-protruding full-length live hinge and extruded tread for superior traction that distinguishes it from its competition. Finished with durable, ergonomically designed handles, this ramp was redesigned to be simply the best in the industry. Top Lip Extension (TLE) compatible.

Available in 2' – 8' lengths with an 800-pound weight capacity.



Features a no-pinch, non-protruding live hinge for safe, silent, and maintenance free operation




Ergonomically designed flexible, non-breakable handles offer convenient carrying and a comfortable grip



Bottom transition plates self-adjust to accommodate uneven terrain

TOP LIP EXTENSION (TLE)

The TLE will increase the reach of the SUITCASE® Signature Series® by extending the top lip from 3" to 9"! The TLE is unique because it can be added or removed as needed, eliminating the need for a large ramp inventory. Shipped as a pair.



Patented



With the use of a TLE, the rear bumper is easily cleared

SUITCASE® Signature Series®

| SPECIFICATIONS | SUITCASE SS2 | SUITCASE SS3 | SUITCASE SS4 | SUITCASE SS5 | SUITCASE SS6 | SUITCASE SS7 | SUITCASE SS8 | TLE |
|----------------------------|---------------|---------------|------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|------------------|
| Made of aluminum | • | • | • | • | • | • | • | • |
| Non-skid surface | • | • | • | • | • | • | • | • |
| Usable size | 2' x 30" | 3' x 30" | 4' x 30" | 5' x 30" | 6' x 30" | 7' x 30" | 8' x 30" | 8½" x 13⅞" x 2⅝" |
| Folded size (L x W x H) | 2' x 15" x 4" | 3' x 15" x 4" | 4' x 15" x 4" | 5' x 15" x 4" | 6' x 15" x 4" | 7' x 15" x 4" | 8' x 15" x 4" | n/a |
| Built-in carrying handle | • | • | • | • | • | • | • | n/a |
| Weight (lbs.) * | 9 | 14 | 18 | 25 | 30 | 36 | 42 | 5 per set |
| Weight capacity (lbs.) | 800 | 800 | 800 | 800 | 800 | 800 | 800 | **500 |
| For wheelchairs & scooters | • | • | • | • | • | • | • | • |
| Usage guidelines | Curbs, 1 step | Curbs, 1 step | Curbs, 1-2 steps | Curbs, 1-2 steps, some minivans | 1-3 steps, some minivans & SUVs | 1-3 steps, some minivans & SUVs | 2-3 steps, some minivans & SUVs | n/a |

* Weights are approximate and may vary. Consult your equipment's owner guide for proper degree of incline. Never exceed its recommendations.

** Ramp weight capacity is reduced when used with the TLE.

