



*Above:  
AL100 lift with swingaway option*

**Easy as 1, 2, 3!**

Harmar's Universal Outside Lifts are one of the few available for vehicles with Class I, II or III hitches. Their quality features, construction and durability make them the most popular outside lifts on the market. The Swing Away option provides easy access to the trunk or lift gate of any vehicle.

**H**armar designs, engineers and produces American made lift, ramp and home access products that enhance people's mobility, independence and quality of life. Our products include a complete line of Auto Lifts, Access Ramps and Home Access Products.

Harmar products are available exclusively through a nationwide network of authorized dealers. Customers receive care from the nation's top mobility dealers. Harmar's dealers will assist in determining the best product to fit your needs and budget. It is imperative to have an authorized Harmar dealer demonstrate, install and service your Harmar lift.

Harmar's lifts, ramps and home access products are designed and built with the "Harmar difference" - they are simpler, stronger, lighter and easier to use.



NEMDA and ADED Member  
We Proudly Serve Veterans



**STEP 1**

Easy to lower with one switch operation.



**STEP 2**

Drive on and secure (automatic on some models).



**STEP 3**

Raise and go!

**Model AL100**

UNIVERSAL SCOOTER LIFT



- America's #1 selling scooter lift
- Fits virtually all scooters
- Automatic hold-down arm
- No scooter modifications
- Large aluminum deck

**Model AL300**

FUSION LIFT



- Carries both power chairs or scooters
- Self tensioning retractors
- No chair or scooter modifications
- Large aluminum deck
- Swing-Away option available

**Model AL300F**

FUSION LIFT



- Carries both power chairs or scooters
- Self tensioning retractors
- No chair or scooter modifications
- Large aluminum deck
- Swing-Away option available



**Model AL300RV**

FUSION LIFT



- Carries both scooters and power chairs
- Self-tensioning retractors
- No chair or scooter modifications
- Manual crank backup
- Swing Away option available

**Model AL301XL**

XL FUSION LIFT



- Carries both larger scooters and power chairs
- Self-tensioning retractors
- No chair or scooter modifications
- Manual crank backup
- Swing Away option available



**Model AL500**

UNIVERSAL POWER CHAIR LIFT

- America's #1 selling power chair lift!
- Fits virtually all power chairs
- Self tensioning retractors
- No chair modifications
- Large aluminum deck
- Swing Away option available

**Model AL550**

AUTOMATIC POWER CHAIR LIFT



- Fits many popular power chairs
- Simple, secure auto-locking feature
- Chair modification required
- Manual crank backup
- Swing away option available

**Model AL560**

UNIVERSAL POWER CHAIR LIFT



- Padded hold down arm secures chair by pressing down on seat.
- Universal design: Works with a wide range of power chairs.
- Platform wheel cradles adjust to provide a custom fit for each chair.
- Hold-down arm adjusts vertically to match your seat height.



**Model AL105**

SWING AWAY OPTION

- Easy access to storage
- Double Safety lock feature
- Allows entry to vehicles with liftgates, rear doors and tailgates

Product	AL100	AL300	AL300F	AL300RV	AL301XL	AL500	AL550	AL560
Platform Size	27.5"x39"-46"	27.5"x39"-46"	27.5"x39"-46"	27.5"x39"-46"	36"x39"-53"	27.5"x39"	27.5"x39"	28.5"x38.5"
Weight Cap.	350lb	350lb	350lb	350lb	350lb	350lb	350lb	350lb
Ins. Weight	95lb	105lb	105lb	110lb	110lb	98lb	112lb	105lb
Works With	Scooters	Both	Both	Both	Both	Power Chairs	Power Chairs	Power Chairs
Hitch Class	II or III	II or III	II or III	III Only	III Only	II or III	II or III	II or III