

## AL500 Universal Power Chair Lift



### FEATURES

- Carries virtually all power chairs
- Four self-tensioning Q'Straint retractors
- No power chair modifications
- Manual crank back up & license plate mount included
- Swing-away option available

### SPECIFICATIONS

Lifting capacity:	350 lbs ( 160 kgs)
Aluminum deck size:	27.5" X 39"
Installed weight:	98 lbs (44.5 kgs)
Hitch height (range):	5"-19"
Hitch Class:	II or III
Warranty:	3-year transferable

### DESCRIPTION

Harmar's durable, American-made AL500 Universal Power Chair Lift is the world's most popular outside platform lift.

The AL500 is easy to operate and built to last for years of trouble-free use. The universal design allows you to safely transport your power chair without making any chair modifications!