

Scout LT 4-Wheel Scooter

HCPCS:



FIGURE A




FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G




FIGURE H

The Scout 4 Wheel comes standard with 12 AH batteries to offer up to 9 miles on a full charge. The Scout is a great travel scooter for those wanting a very lightweight scooter that can go the distance for a full day of activities. Your needs are sure to be met with all of the standard features on the Scout: a Delta Tiller for one-handed operation, adjustable armrests, a generously-sized front basket for storage of personal items, anti-tip wheels, and low-maintenance, flat-free tires.

- A fun, easy to disassemble 4-wheel scooter for the budget conscious shopper.
- Includes easily interchangeable Red and Blue shrouds
- For SCOUTLT4, Up to 9 miles per charge
- For SCOUTLT4-EXT, Up to 14 miles per charge
- Up to 9 miles per charge
- Non-marking, solid flat-free tires
- Padded Delta tiller

PRODUCTS	Description
SCOUTLT4	Scout LT 4-Wheel Scooter, Red/Blue, 1/ea
SCOUTLT4-EXT	Scout LT 4-Wheel 20AH Battery, Red/Blue
SPECIFICATIONS	Description
Warranty	24 Month Limited Warranty on Electronics
Product Weight	100 lb
Product Dimensions	42.3" (L) x 19.7" (W) x 36.2" (H)
Seat Dimensions	17.3" (D) x 14" (W)
Weight Capacity	300 lb
Top Speed	4 mph
Battery Charger	2A Offboard
Climbing Angle	6°
Color	Red & Blue
Controller	PG 45A
Ground Clearance	3.3"
HCPCS	K0800
Motor	24V x 250W x 170
Seat Type	Stadium
Turning Radius	53"
Warranty on Batteries	12 Month Warranty on Batteries
Warranty on Frame	Limited Lifetime Warranty on Frame
Wheel Size	Front: 8" x 2" Rear: 8" x 2"

SPECIFICATIONS

Description

CA Prop 65

[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.]