

Triaxe Sport Specifications



Model	Triaxe Sport T3045
Overall Length	34"
Overall Width	20 1/2 "
Tires	8" Solid Front Tire 9.5 Pneumatic Rear Tires
Maximum Speed	Up to 9mph
Safe Gradient	0-18 degrees
Range	Up to 35 miles
Turning Radius	33"
Weight (without battery pack)	49lbs
Color	Metallic Red, Light Metallic Blue
Brake	Hand brake controlling rear disc brakes
Drive System	Front Wheel Drive, Brushless Motor
Handlebar	Steel with Rubber Hand Grips
Maximum Weight Capacity	Up to 350lbs
Seat Width	13" standard, 18" and 20" optional
Seat Height	20.5"
Battery	48V 20amp – 8lbs
Motor	56V -350W
Battery Charger	48V – 2AH

Revised 4/16