

There is Power Mobility. And then there's

POWERFUL MOBILITY



transport
AF

Call 1-800-944-1997 or visit us at www.evrider.com

There is Power Mobility. And then there's

POWERFUL MOBILITY



WEIGHT CAPACITY	250 lbs
DRIVE WHEEL	8" x 2", Solid Tire, Flat Free
FRONT WHEEL	7" x 1.6", Solid Tire, Flat Free
REAR CASTER (ANTI-TIP-PER)	yes, standard
MAX SPEED	up to 3.8 mph *
BATTERY SPECIFICATIONS	24V 8.7Ah Lithium Battery x1 (included)
MAXIMUM TRAVEL RANGE	up to 6.8 miles*

CHARGER TYPE	2 Amp 24 VDC, input: 110/240 Vac, Off Board Li-Ion Charger
CONTROLLER TYPE	Dynamic R-Series 50A
MOTOR TYPE	180 watt 24 volt
LENGTH	33.8"
WIDTH	16.5"
OVERALL HEIGHT	32.8"
WEIGHT	44 lbs
SUSPENSION	none

TURNING RADIUS	31"
SEAT HEIGHT	20"
WHEELBASE	25"
GROUND CLEARANCE	1.4"
OVERALL DIMENSIONS	33.8" x 16.5" x 32.8"
FOLDED DIMENSIONS	28.7" x 17.3" x 16.5"

(*) SPECIFICATION PERFORMANCE IS DETERMINED IN THE LABORATORY SETTING BASED ON IDEAL CONDITIONS (MINIMAL FRICTION, NO WIND AND CONTROLLED TEMPERATURE AND PRESSURE WITH AN AVERAGE STANDARD WEIGHT). ACTUAL PERFORMANCE WILL BE DIFFERENT THAN RESULTS OBTAINED DURING LABORATORY TESTING.

Call 1-800-944-1997 or visit us at www.evrider.com