

**Avenger** 

4 wheel

Seat Color Charcoal Vinyl

Body Colors Silver Re

> The Golden Avenger four-wheel, heavy duty scooter provides an impressive blend of rugged outdoor maneuverability while giving a smooth, comfortable ride over a variety of

nature's most challenging terrains. Featuring a 500 lb. weight capacity, the Avenger comes standard with NF22 batteries providing a range of 18 miles.





GA-541

## Avenger 4 wheel **Specifications**

*	
Model	GA541
Medicare Code	Not Classified
Weight Capacity	500 lbs
Drive Wheels	Rear
Maximum Speed	7 mph
Operating Range 1	18 miles
Ground Clearance to Center D	
Ground Clearance to Transaxle	
Turning Radius	67"
Type Batteries	2-NF22
Colors	Silver, Red
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	PG 110 Amp
Dimensions:	
Length	57.5"
Width	24.5"
Height (ground to top	
of back of seat)	42.5" - 45.5"
Height (ground to top	
of back of headrest)	50" - 56"
Ground to Top of Seat	25" - 28"
Ground to Top of Deck	6.5"
	18.75" - 21.75"
Number of Seat Height Posit	
Size of Increments	1"
Front Axle to Rear Axle	38.5"
Weight of Unit: (assembled	-
Front Half	68 lbs
Rear Half	68 lbs
Base	N/A
Seat w/Arms	55 lbs
Batteries	37 lbs ea
Tires:	Pneumatic
Front	13" 13"
Rear	
Casters	N/A
Rear Anti-Tippers	3"
Standard Seat: Slide Seat	Captain's Seat N/A
	21" - 24"
Back Height - no headrest Back Height with headrest	27.5" - 31.5"
•	27.5" - 31.5" 20" X 17"
Width x Depth Color	
	Gray Vinyl
Warranty:	:fat:=
Frame Drive Train	Lifetime
Drive Train	2 Years
Electronics <sup>2</sup>	2 Years

<sup>&</sup>lt;sup>1.</sup> Battery range will vary due to rider weight, drive surface and terrain. Foam-filled tires are optional. The operational mileage will be less if scooter is equipped



Easy to

Battery Meter

Easy to

Charger

Reach

Plug

Read

## **Features**

- 500 lb weight capacity
- Speeds up to 7 mph
- Bright safety lights
- 22 NF batteries
- Up to 18 Mile Range
- 4-Way Flashers Great Warranty

- Removable basket









Autolight

Key Switch

Rubber

Grips

with foam filled tires.
<sup>2</sup> Electronics warranty excludes batteries.