

FULL SIZE PORTABLE SCOOTER



U1 Battery for Maximum Range

Seat Color
Black Vinyl

Body Colors
Blue
Red



3 wheel GB118
4 wheel GB148

The all new Golden BuzzAround Extreme Scooter features the portability and easy disassembly of a lightweight scooter with the amazing benefits of a full size scooter! The BuzzAround Extreme Scooters come standard with 2-U1 batteries - industry leading and the largest battery in this category - for maximum operating range up to 18 miles. The Extreme also features a brand new ultrabright LED headlight and taillight, 330 lb. weight capacity, larger 9" tires, and comfort suspension springs for the smoothest ride. Go the extra mile on a BuzzAround Extreme Scooter!





FULL SIZE PORTABLE SCOOTER

Wireless disassembly in 4 easy steps!

Item	GB118	GB148
Description	EX	EX
Medicare Code	Not Coded	Not Coded
Weight Capacity	330 lbs.	330 lbs.
Drive Wheels	Rear	Rear
Maximum Speed	5 mph	5 mph
Operating Range up to ¹	18 miles	18 miles
Under Deck Clearance	4"	4"
Turning Radius	42"	56"
Type Batteries	2-U1	2-U1
Colors	Red, Blue	Red, Blue
Freewheel Mode	Yes	Yes
Electronic Speed Control	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes
Charger	Off-board	Off-board
Controller	90 Amp Dynamic R-Series	90 Amp Dynamic R-Series

Dimensions		
Length	45"	44"
Width	22"	22"
Height (ground to top of back of seat)	35"-37"	35"-37"
Ground to Top of Seat	22.75"-24.75"	22.75"-24.75"
Ground to Top of Deck	6.75"	6.75"
Deck to Top of Seat	16"-18"	16"-18"
Number of Seat Height Positions	3	3
Size of Increments	1"	1"
Front Axle to Rear Axle	34.7"	33.8"

Weight of Unit: (assembled)		
Front Half	148 lbs.	155 lbs.
Rear Half	46 lbs.	53 lbs.
Seat w/Arms	30 lbs.	30 lbs.
Batteries	23 lbs.	23 lbs.
	24 lbs. per pack	24 lbs. per pack
Tires:		
Front	Solid	Solid
Rear	9"	9"
Rear Anti-Tip Wheels	9"	9"
	2"	2"

Standard Seat		
Seat Back Height - no headrest	Stadium Style	Stadium Style
Width x Depth	16"	16"
	18" x 16"	18" x 16"
Color	Black Vinyl	Black Vinyl
Optional Seat:		
Back Height - no headrest	Stadium Style	Stadium Style
Width x Depth	16"	16"
	20" x 18"	20" x 18"
Color	Black Vinyl	Black Vinyl

Warranty:		
Frame	Lifetime Warranty	Lifetime Warranty
Drive Train	2 Year Warranty	2 Year Warranty
Electronics ²	2 Year Warranty	2 Year Warranty
Batteries	1 Year Warranty	1 Year Warranty

¹ Battery range tested at 200 lbs. and will vary due to rider weight, drive surface and terrain.

² Electronics warranty excludes batteries.



Step 1: Pop off the seat!



Step 2: Pull out the battery pack!



Step 3: Fold down tiller and lock in place!



Step 4: Pull forward on the drivetrain release lever to separate the front and rear frame sections!



Disassembles into 6 pieces!

4 wheel model shown

Features

- 330 lb. weight capacity
- 4" ground clearance
- Full size & portable
- Larger 90 amp controller
- Larger 320 watt motor
- Larger 9" tires
- U1 battery for maximum range
- Ultrabright LED tail light
- Comfort suspension springs
- Full size, yet portable & easy to disassemble



A battery charging port is conveniently located on the tiller under the control panel.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

