



Golden Scooters

SUPERIOR DESIGNS | FIRST CLASS LUXURY

SUMMER 2016



 **GOLDEN**

BUZZAROUND XLTM SERIES

COMPACT PORTABLE

The next generation of the popular Buzzaround scooter series is now available! The Buzzaround XL offers more standard features to make your riding experience better than ever. Featuring a wrap around Delta Tiller, a new LED battery gauge, and a new ultra bright LED angle-adjustable headlight, the new Buzzaround XL is also longer to offer more foot room, with a great turning radius. The Buzzaround XL disassembles in seconds with just 4 easy steps making transporting it from one destination to another a breeze!





GB106D

BUZZaround®
lite

BUZZAROUND LITE

- 250 lb. weight capacity
- Delta Tiller standard
- Disassembles into 5 lightweight pieces
- Top speed 4 mph
- Three sets of interchangeable shroud panels **INCLUDED FREE:** Red, Blue and Sahara

***"Life is just better
on a Buzzaround!"***

(Buzz XL model shown below)





GB117D



GB147D

BUZZAROUND XL MODELS FEATURE:

- 300 lb. weight capacity
- Ultra Bright LED headlight standard
- Charging port below the control panel
- 3 inch ground clearance
- Red and Blue Shroud panels included free
- Disassembles quickly and easily for transport

WARRANTIES:

- FRAME:** Lifetime Warranty
- DRIVE TRAIN:** 2 Year Warranty
- ELECTRONICS:** 2 Year Warranty
- BATTERIES:** 1 Year Warranty



BUZZAROUND XL HD

GB117H



BUZZAROUND XL HD

GB147H

ALL BUZZAROUND XLHD MODELS FEATURE:

- Increased weight capacities up to 325 lbs.
- Ultra Bright LED headlight standard
- Charging port below the control panel
- 3 inch ground clearance
- Red and Blue Shroud panels included free
- Disassembles quickly and easily for transport
- 22Amp batteries standard
- Operating range up to 14 miles per charge

WARRANTIES:

FRAME: Lifetime Warranty
DRIVE TRAIN: 2 Year Warranty
ELECTRONICS: 2 Year Warranty
BATTERIES: 1 Year Warranty

**BUZZY'S
AROUND XLS™**
Comfort Spring Suspension



GB117S

**BUZZY'S
AROUND XLS™**
Comfort Spring Suspension



GB147S

BUZZAROUND XLS MODELS FEATURE:

- 300 lb. weight capacity
- Ultra Bright LED headlight standard
- Charging port below the control panel
- 3 inch ground clearance
- Red and Blue Shroud panels included free
- Disassembles quickly and easily for transport
- NEW Comfort Spring Suspension for a smoother ride

WARRANTIES:

FRAME: Lifetime Warranty
DRIVE TRAIN: 2 Year Warranty
ELECTRONICS: 2 Year Warranty
BATTERIES: 1 Year Warranty

BUZZAROUND XLS[™] HD



GB117Z



BUZZAROUND XLS[™] HD

GB147Z

**BUZZAROUND XLSHD
MODELS FEATURE:**

- Increased weight capacity up to 325 lbs.
- Ultra Bright LED headlight standard
- Charging port below the control panel
- 3 inch ground clearance
- Red and Blue Shroud panels included free
- Disassembles quickly and easily for transport
- 22 Amp batteries standard
- Operating range up to 14 miles per charge
- NEW Comfort Spring Suspension for a smoother ride

WARRANTIES:

- FRAME:** Lifetime Warranty
- DRIVE TRAIN:** 2 Year Warranty
- ELECTRONICS:** 2 Year Warranty
- BATTERIES:** 1 Year Warranty

BUZZ[™]AROUND EX[™]

EXTREME

FULL-SIZE PORTABLE

The all new Golden Buzzaround Extreme Scooter features the portability and easy disassembly of a lightweight scooter with the amazing benefits of a full size scooter. Both the three and four-wheel BuzzAround Extreme Scooters come standard with 2-U1 batteries – industry leading and the largest battery in this category – for maximum operating range up to 18 miles. The Extreme also features a brand new ultrabright LED headlight and taillight, and comfort suspension springs for the smoothest ride.

Go the extra mile on a Buzzaround Extreme Scooter!

BUZZAROUND EXTREME FEATURES:

- 330 lb. weight capacity
- 4" ground clearance
- Full size & portable, yet easy to disassemble
- Larger 90 amp controller
- Larger 320 watt motor
- Larger 9" tires
- 2-U1 batteries for maximum range
- Ultrabright LED headlight and tail light
- Comfort suspension springs



WARRANTIES:

FRAME: Lifetime Warranty

DRIVE TRAIN: 2 Year Warranty

ELECTRONICS: 2 Year Warranty

BATTERIES: 1 Year Warranty



GB118



GB148



LiteRider^{GOLDEN}[®]

LARGE PORTABLE

The unique tubular design of the LiteRider scooters gets you wherever you need to go in a style all their own. The only scooters to come with two cargo baskets standard, the LiteRiders are easy to disassemble, and maneuver smoothly indoors and out.

LiteRider^{GOLDEN}[®]
LiteRider GL111D



LiteRider^{GOLDEN}[®]
LiteRider GL141D

THE GL111D AND GL141D FEATURES:

- 300 lb. weight capacity
- Wraparound Delta Tiller standard
- Adjustable tiller standard
- Easy to read control panel with backlit battery gauge
- Two storage baskets standard
- Comfortable 17"x16" stadium style seat
- Adjustable armrests
- 2-22 Ah Batteries; increased operating range up to 12 miles
- Completely wireless disassembly
- 9" front and rear tires
- New Ultrabright LED headlight and rugged plastic basket standard!

LiteRider^{GOLDEN}
Mini

LiteRider GL109D



- 250 lb. weight capacity
- New seat design
- Operating range up to 8 miles
- Turning radius: 36.5"
- 2-12 Ah Batteries
- 8" front and rear tires
- Two storage baskets standard
- Adjustable armrests



LiteRider^{GOLDEN}
ENVY

LiteRider GP162

Golden's RETAIL CASH POWER WHEELCHAIR features:

- 300 lb. weight capacity
- 2-22 Ah Batteries; increased operating range up to 15.5 miles
- Turning radius: 26"
- Single post seat mount
- Disassembles easily for transport



GOLDEN COMPANION

FULL-SIZE LUXURY

The Golden Companion scooters are your ticket to ride first class every day! Allow the luxurious, full-size Companions to enhance your life and take you to the places you've always had trouble getting to. Shop in the mall, attend sporting events and enjoy the great outdoors with the unmatched comfort of your favorite new Companion! All Companion models feature:

- All new design for easy disassembly
- Wraparound Delta Tiller standard
- Adjustable LED headlight
- EZ Reach tiller adjustment lever
- Auto lock tiller
- Easy to read control panel
- Easy access charging port
- Backlit battery gauge
- Black, non-marking, low profile tires on sturdy steel rims
- Comfortable, fully adjustable padded armrests
- Under-carpet troubleshooting guide
- Removable shopping basket with handle
- Four mobility transportation tie downs



Companion GC240

Companion GC440

Companion GC340

golden **Avenger**



Avenger GR541D

When you want a rugged, heavy duty scooters to tackle the challenging terrain of the great outdoors, look no further than the Golden Patriot™ and Avenger™. Both feature oversized tires with sporty mag wheels, a full-size captain's seat for all day riding comfort, and an impressive operating range of up to 18 miles.

The Avenger model features:

- 500 lb. weight capacity
- Adjustable tiller for riding comfort
- Deluxe, reclining high back captain's seat
- Large, 13" front and rear tires
- Adjustable front headlight and four-way flashers
- Large front cargo basket standard
- 18-mile operating range

Golden **Patriot**™

The Patriot™ model features:

- Full front and rear lighting package with tail lights, headlights and flashers
- Four-wheel independent suspension for smoothly traversing outdoor terrain
- Rugged rubber mat and anti-skid material on the floor deck
- Convenient water bottle holders for two cold beverage bottles
- Full control panel featuring safety horn, flashers, lights, turn signals, speed knob and backlit battery gauge
- Black, low profile, non-marking solid tires with chrome hubcaps



Patriot GR575D

AVAILABLE ACCESSORIES

**Only available on Companion models.*



Oxygen Holder



Crutch Holder



Rear Basket



Safety Flag



Comfort Ride*



Walker Holder



Power Elevating Seat*



Cup Holder



Quad Cane Holder



Pack'n'Go



Cell Phone Holder

SPECIFICATIONS*

| Item | Description | Length | Width | Front Tire | Rear Tire | Weight Capacity (lbs.) | Top Speed (mph) | Under Deck Clearance | Operating Range** (miles) | Turning Radius | Batteries (2) | Weight w/ Batteries | Front Half Weight | Rear Half Weight | Seat Weight w/ Arms | Handlebar Piece | Freeheel Mode | Electronic Speed Control | Charger |
|--------|----------------------|--------|--------|------------|-----------|------------------------|-----------------|----------------------|---------------------------|----------------|---------------|---------------------|-------------------|------------------|---------------------|-----------------|---------------|--------------------------|-----------|
| GB106D | Buzz Lite (3-wheel) | 36.5" | 21.5" | 8" | 8" | 250 | 4 | 1.75" | 8 | 31" | 2-12Ah | 96 | 29 | 26 | 22 | 29 | Yes | Yes | Off-board |
| GB117D | Buzz XL (3-wheel) | 42.5" | 21.5" | 8" | 8" | 300 | 4 | 3" | 8' | 37" | 2-12Ah | 101 | 35 | 25 | 20 | 35 | Yes | Yes | Off-board |
| GB117S | Buzz XLS (3-wheel) | 42.5" | 21.5" | 8" | 8" | 300 | 4 | 3" | 8' | 37" | 2-12Ah | 103 | 35 | 28 | 21 | 35 | Yes | Yes | Off-board |
| GB117H | Buzz XLHD (3-wheel) | 42.5" | 21.5" | 8" | 8" | 325 | 4 | 3" | 14' | 37" | 2-22Ah | 110 | 35 | 25 | 21 | 35 | Yes | Yes | Off-board |
| GB117Z | Buzz XLSHD (3-wheel) | 42.5" | 21.5" | 8" | 8" | 325 | 4 | 3" | 14' | 37" | 2-22Ah | 114 | 35 | 28 | 22 | 35 | Yes | Yes | Off-board |
| GB147D | Buzz XL (4-wheel) | 41.5" | 21.5" | 8" | 8" | 300 | 4 | 3" | 8' | 47" | 2-12Ah | 109 | 43 | 25 | 20 | 43 | Yes | Yes | Off-board |
| GB147S | Buzz XLS (4-wheel) | 41.5" | 21.5" | 8" | 8" | 300 | 4 | 3" | 8' | 47" | 2-12Ah | 111 | 43 | 28 | 21 | 43 | Yes | Yes | Off-board |
| GB147H | Buzz XLHD (4-wheel) | 41.5" | 21.5" | 8" | 8" | 325 | 4 | 3" | 14' | 47" | 2-22Ah | 118 | 43 | 25 | 21 | 43 | Yes | Yes | Off-board |
| GB147Z | Buzz XLSHD (4-wheel) | 41.5" | 21.5" | 8" | 8" | 325 | 4 | 3" | 14' | 47" | 2-22Ah | 122 | 43 | 28 | 22 | 43 | Yes | Yes | Off-board |
| GB118 | Buzz XLEX (3-wheel) | 45" | 22" | 9" | 9" | 330 | 5 | 4" | 18' | 42" | 2-U1 | 148 | 46 | 30 | 23 | 46 | Yes | Yes | Off-board |
| GB148 | Buzz XLEX (4-wheel) | 44" | 22" | 9" | 9" | 330 | 5 | 4" | 18' | 56" | 2-U1 | 155 | 53 | 30 | 23 | 53 | Yes | Yes | Off-board |
| GL111D | LiteRider (3-wheel) | 41.5" | 22" | 9" | 9" | 300 | 5 | 3.75" | 12 | 39" | 2-22Ah | 119 | 44 | 28 | 22 | 44 | Yes | Yes | Off-board |
| GL141D | LiteRider (4-wheel) | 43.25" | 22" | 9" | 9" | 300 | 5 | 3.75" | 12 | 48" | 2-22Ah | 127 | 52 | 28 | 22 | 52 | Yes | Yes | Off-board |
| GL109D | LiteRider Mini | 37.75" | 19.75" | 8" | 8" | 250 | 4 | 2.75" | 8.4' | 36.5" | 2-12Ah | 101 | 37 | 22 | 24 | 37 | Yes | Yes | Off-board |
| GP162 | LiteRider Envy | 33" | 22.5" | N/A | 9" | 300 | 3.5 | 2.5" | 15.5' | 26" | 2-22Ah | 120 | 34 | 35 | 22 | 35 | Yes | Yes | Off-board |
| GC240D | Companion (3-wheel) | 42.5" | 23.5" | 9" | 9" | 300 | 4.5 | 3.75" | 15.5 | 39" | 2-U1 | 166 | 46 | 48 | 28 | 48 | Yes | Yes | Off-board |
| GC340D | Companion (3-wheel) | 45.5" | 23.5" | 9" | 10" | 350 | 4.5 | 4.5" | 16.5 | 40.5" | 2-U1 | 168 | 48 | 48 | 28 | 48 | Yes | Yes | Off-board |
| GC440D | Companion (4-wheel) | 47.5" | 23.5" | 10" | 10" | 350 | 4.5 | 4.5" | 15 | 53.5" | 2-U1 | 183 | 63 | 48 | 28 | 63 | Yes | Yes | Off-board |
| GA541D | Avenger | 57.5" | 24.5" | 13" | 13" | 500 | 7 | 4" | 18 | 67" | 2-NF22 | 265 | 68 | 68 | 55 | 68 | Yes | Yes | Off-board |
| GR575D | Patriot | 57" | 28" | 13" | 13" | 400 | 8 | 3.5" | 18 | 76.5" | 2-75A h | 314 | Base 158 | N/A | 50 | 160 | Yes | Yes | Off-board |

* ALL WEIGHT SHOWN IN LBS.

**Battery range will vary due to rider weight, drive surface and drive slope. Range was tested at maximum weight.

1. Range tested at 200 lbs.

Brakes: All Golden Technologies scooters come standard with Dynamic/Electromechanical braking systems. Rear wheel drive: All Golden Technologies scooters come standard with a sealed transaxle, 24 volt motor.

Choose from a Complete Line of Lifestyle Solution Products by Golden



POWER WHEEL CHAIRS

- Compass
- Compass HD
- LiteRider Envy
- LiteRider Mini
- Alanté
- Alanté Sport



POWERED LIFT AND RECLINE CHAIRS

- Value Series
- Traditional Series
- Signature Series
- Comforter Series
- Comforter Wide Series
- MaxiComforter Series

Visit Your Local Authorized Golden Dealer:



www.goldentech.com