



5th Generation Companion Scooters & ALL-NEW Companion HD

Product Overview



Golden's Companion Scooter is Tried-&-True!



Golden's Popular Companion Scooter Has a Sporty New Look & Upgraded Comfort Features!

The **Next-Generation Companion™** combines a modern design with high performance features to offer the ultimate indoor or outdoor driving experience! Built for stability and maneuverability, it offers the most leg and foot room in its class with upgraded comfort features for a dependable full-time, yet portable scooter.



GOLDEN **Companion™**

GC240

Mid-Size 3-Wheel Scooter

MAP \$2,156

MSRP \$2,756

Weight Capacity: 350 lbs.

Weight of Unit (assembled): 170.8 lbs.

NEW FOR 2024!

- **NEW High Back Stadium Seat with Quilted Diamond Stitch**
- **NEW Azure Blue & Crimson Red Color Shrouds**



 **GOLDEN**

GOLDEN Companion™



GC340 – 5th Generation!

Full-Size 3-Wheel Scooter

MAP \$2,486

MSRP \$2,864

Weight Capacity: 400 lbs.

Weight of Unit (assembled): 207 lbs.

- **NEW Tiller Design with Digital Dash & Console**
- **NEW Tiller Mounted USB Charging Port**
- **NEW High Back Stadium Seat with Quilted Diamond Stitch**
- **NEW Full Front & Rear Suspension**
- **NEW Larger Optional Battery Size for Increased Range**
- **NEW Galactic Grey, Azure Blue & Crimson Red Color Shrouds**
- **3.5 Amp Super Charger for Maximum Battery Life**



GOLDEN **Companion™**

GC440 – 5th Generation!

Full-Size 4-Wheel Scooter

MAP \$2,589

MSRP \$3,135

Weight Capacity: 400 lbs.

Weight of Unit (assembled): 222 lbs.

- **NEW Tiller Design with Digital Dash & Console**
- **NEW Tiller Mounted USB Charging Port**
- **NEW High Back Stadium Seat with Quilted Diamond Stitch**
- **NEW Full Front & Rear Suspension**
- **NEW Larger Optional Battery Size for Increased Range**
- **NEW Galactic Grey, Azure Blue & Crimson Red Color Shrouds**
- **3.5 Amp Super Charger for Maximum Battery Life**



GOLDEN **Companion™**

NEW GC340 & GC440 Features



Luxurious High Back Stadium Seat
with Quilted Diamond Stitch



Digital Dash with Percentage Battery Life Indicator



Full Front & Rear Suspension



2 Convenient Drink Holders



Sporty Black Wheels on
Galactic Grey Models

Introducing the GC540 Companion HD Scooter!

GOLDEN
Companion[™]
HD

The **ALL-NEW Golden Companion HD** 3-wheel scooter offers a modern design with rugged construction built for stability and maneuverability with front anti-tip wheels and a 51" turning radius. It offers a weight capacity of 500 lbs., rear suspension for a smooth ride, and a powerful 120 Amp controller for outstanding climbing ability. It features front and rear LED lighting, an adjustable tiller with auto-lock, and a luxury high-back Captain's seat for maximum comfort when on the go, plus durable vibrant red, blue or grey color shrouds for a sporty look. With solid, non-marking tires, an operating range of 16 miles on 2-NF22 batteries, and a maximum speed of 6 mph, the Companion HD offers everything you need to enjoy the great outdoors!



 **GOLDEN**

GOLDEN CompanionTM HD



GC540 – ALL-NEW!

Heavy-Duty 3-Wheel Scooter

MAP \$2,999

MSRP \$4,799

Weight Capacity: 500 lbs.

Weight of Unit (assembled): 225.5 lbs.

- 500 lb. Weight Capacity
- Extra Wide High Back Seat
- Tiller Mounted USB Charging Port
- Rear Suspension
- Galactic Grey, Azure Blue & Crimson Red Color Shrouds
- 8 Amp Super Charger for Maximum Battery Life

Medicare Code: K0808



GOLDEN Companion™ HD

NEW GC540 Features



Digital Dash with Percentage
Battery Life Indicator



Rear Suspension for a Smooth Ride



2 Convenient Drink Holders



Front & Rear LED Lights & Side
Reflectors for Safety





Marketing Assets:

- Banner
- Ad Slicks
- Sell Sheets
- Product Images
- Social Media Toolkit

**Order Your New
Companion
Scooters Today!**

REVAMPED & READY!

**GOLDEN
Companion™**



AVAILABLE IN MID-SIZE, FULL-SIZE
& NEW HEAVY-DUTY OPTIONS!

The Next-Generation Companion™
is a sleek, high-performance scooter designed for both
indoor and outdoor adventures! With top-notch stability,
maneuverability, and spacious leg and foot room,
it ensures a reliable full-time experience.

New Features:

- Sporty Tiller Design with Two Convenient Drink Holders
- Digital Dash with Percentage Battery Life Indicator
- Water-Resistant Buttons on the Control Panel
- Luxurious High Back Stadium Seat with Quilted Diamond Stitch
- Larger Optional Battery Size for Increased Range (GC340 & GC440)
- Full Front & Rear Suspension on GC340 & GC440 (Rear Suspension on GC540)
- 3.5 Amp Super Charger for Maximum Battery Life on GC340 & GC440 (8 Amp on GC540)



GOLDEN