

S549--ARMU/ S549--SBMU


| | |
|---|----------------------------------------------|
| ① | Seat ass'y |
| | 32200239 (SEAT W/ ARMRESTS,ASS'Y) |
| | 43701035 (SEAT,CAPTAIN,SMALL) |
| | 31005013 (ARMREST ASS'Y,LH,BLACK) |
| | 31005014 (ARMREST ASS'Y,RH,BLACK) |
| | 43704001 (ARMREST PAD,DESK ARM,BLACK) |
| ② | Battery Pack Ass'y. |
| | 33700174 (BATTERY PACK,15AH) |
| ③ | Rear Shroud Ass'y. |
| | 30510927 (SHROUD ASS'Y,REAR,CANDY APPLE RED) |
| | 30515705 (SHROUD ASS'Y,REAR, METALLIC BLUE) |
| ④ | Drive Wheel Ass'y |
| | 32600422 (TIRE,PU,MAG WHEEL,DRIVE WHEEL) |
| ⑤ | Motor & Transmission. |
| | 62000054 (TRANSAXLE ASS'Y/ M630CA-4) |
| | 62100092 (MOTOR ASS'Y/ M6CA4) |
| | 62200129 (TRANSAXLE/M630) |

| | |
|----|-----------------------------------------------|
| ⑥ | Controller |
| | 64520005(CONTROLLER,P&G S-DRIVE45) |
| ⑦ | Tiller Ass'y. |
| | 33510852 (TILLER,ASS'Y, CANDY APPLE RED) |
| | 33515832 (TILLER,ASS'Y, METALLIC BLUE) |
| ⑧ | Basket |
| | 41700058 (BASKET) |
| ⑨ | Front Shroud Ass'y. |
| | 30510908 (SHROUD ASS'Y,FRONT,CANDY APPLE RED) |
| | 30515998 (SHROUD ASS'Y,FRONT,METALLIC BLUE) |
| ⑩ | Front Wheel Ass'y. |
| | 33100676 (TIRE,PU,MAG WHEEL,FRONT WHEEL) |
| 11 | Charger unit |
| | 65100001(CHARGER,24V/1.5AHP-1201-2-OFF BOARD) |