

PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

GO GO[®] SPORT

3 & 4-WHEEL SCOOTER

FEATURES,
PERFORMANCE,
AND
TRAVEL READY

The Go-Go[®] Sport delivers high-performance and convenience on the go. With a 325 lb. weight capacity, a charger port in the tiller, standard front LED lighting and a maximum speed up to 4.7 mph., the Go-Go[®] Sport is feature rich and travel ready.



Go-Go[®] Sport 3-wheel



Go-Go[®] Sport 4-wheel



Wraparound delta tiller

pridemobility.com



SPECIFICATIONS

GO-GO® SPORT 3W: SC73 (K0801) 4W: SC74

Weight capacity:	325 lbs.
Maximum speed ¹ :	Up to 4.7 mph
Ground clearance ² :	3" at frame, 1.5" at motor
Turning radius:	3W: 36" 4W: 50.5"
Overall length ² :	3W: 42" 4W: 43"
Overall width ² :	21.25"
Front tires:	3W: 3" x 8" solid 4W: 3" x 9" solid
Rear tires:	3" x 9" solid
Range per charge ^{1,4} :	3W: Up to 12.7 miles 4W: Up to 12 miles
Total weight without batteries ⁵ :	3W: 92.8 lbs. 4W: 99.2 lbs.
Heaviest piece when disassembled:	3W: 35 lbs. (front section) 4W: 44 lbs. (front section)
Standard seating:	Type: Foldable, compact, vinyl-covered molded plastic Weight: 23 lbs. Material: black vinyl Dimensions: Width: 17", Depth: 17"
Motor output:	24-volt DC motor, sealed transaxle
Brakes:	Regenerative and electromechanical
Battery requirements ^{3,5} :	(2) 12V deep-cycle (standard) Size: (2) 18 AH, 14 lbs. each Weight: 29.5 lbs. per battery pack
Battery charger:	Off-board, 2 amp
Color:	2 sets of exciting colors with every scooter!



Red



Blue

- ¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- ² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- ³ AGM or gel cell type required.
- ⁴ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance test conducted at maximum weight capacity.
- ⁵ Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

FEATURES

- Feather-touch disassembly permits simple frame separation with only one hand and easily disassembles into five super lightweight pieces for convenient transport and storage
- Black, non-scuffing tires exclusively from Pride®
- Includes two sets of easily changeable colored panels in red and blue
- Front frame-mounted seat post offers maximum stability
- Microprocessor-based controller offers optimal power management and added safety features
- Standard front tiller-mounted basket included
- Deck basket for additional storage (4-wheel only)
- Convenient off-board dual voltage charger can charge battery pack on-board or off-board
- Charger port in tiller

ACCESSORIES

- Rearview mirror
- Cane/crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Wishbone crutch holder
- Double crutch holder
- Weather cover
- Saddlebag
- Cell phone holder
- Lap belt
- Under seat storage (assembly required)
- Additional 18AH battery box

*Military Patches Available



LIFT FINDER

Mobility Made Easy. Pride offers a variety of lifts to transport your mobility product with ease. To find a compatible lift for your vehicle and mobility product, use our new Lift Finder application at pridemobility.com/lift-finder/