

# SPRINTER™ XL4 DELUXE 889-4XLS

Introducing the **Sprinter XL4 Deluxe**, fully equipped with many features such as full suspension, full light package and delta style tiller. The four-wheel design and suspension lets the **Sprinter XL4 Deluxe** drive with ease. The **Sprinter XL4 Deluxe** allows a steady ride with lots of flair.

## Available Colors



Burgundy



- **500 lb. weight capacity**
- **20" deluxe captain seat that slides, swivels and reclines**
- **2 HP motor**
- **Full light package**
- **Full suspension**
- **12" wheels**
- **Delta style adjustable tiller**
- **Up to 25 mis. range**
- **20 speed settings**
- **Emergency hand brake**
- **Basket included**

Overall Dimensions (L1 x W x H)	in	53 x 27 x 34
Number, Size of Tires (front, rear)	Front (Rear)	2, 12" (2, 12")
Suggested User Weight (on a level road)	lb (kg)	Max. 500 (227)
Weight of Heaviest Part (w/o baskets and seat)	lb (kg)	115 (52)
Total Weight (w/ batteries, basket, seat)	lb (kg)	240 (109)
Each Battery Module Weight	lb (kg)	50 (23)
Battery Module Capacity	_V_Ah x Pcs	12V33-75Ah x 2
Charger		Off Board (8A)
Output Power of Motor	hp (W)	2 (1491)
Maximum Speed <sup>2</sup>	mph (km/hr)	5 (8)
Maximum Climbing Angle / Safe Climbing Angle (based on user weight of 90kg/200lb)	% (°)	14 (8)*
Range (per charge) (After the battery and mechanical moving parts fully break in)	mile (km)	25 (40)
Turning Radius	in (mm)	72 (1829)
Ground Clearance	in (mm)	4 (102)
Handle Bar	Type	Delta

**SHOPRIDER**

Toll-Free: 800.743.0772  
www.shoprider.com