

SPRINTER™ XL4

8889B-4

With a 350 lb. weight capacity, the **Sprinter XL4** is a perfect heavy duty scooter. The **Sprinter XL4** comes with a deluxe captain seat, chrome bumper and 12" wheels for added comfort. The **Sprinter XL4** is highly maneuverable and is loaded with standard features at a value price.

Available Colors



Burgundy



- 350 lb. weight capacity
- Front chrome bumper
- Slide and swivel deluxe captain seat
- Front headlight package
- 12" wheels for added safety
- Adjustable tiller
- Rear light package w/blinkers/hazards avail.

Overall Dimensions (L1 x W x H)	in	52 x 23 x 30-42
Number, Size of Tires (front, rear)	Front (Rear)	2, 12" (2, 12")
Suggested User Weight (on a level road)	lb (kg)	Max. 350 (159)
Weight of Heaviest Part (w/o baskets and seat)	lb (kg)	108 (49)
Total Weight (w/ batteries, basket, seat)	lb (kg)	181 (82)
Each Battery Module Weight	lb (kg)	32 (12)
Battery Module Capacity	_V_Ah x Pcs	12V33-40Ah x 2
Charger		Off Board (3A)
Output Power of Motor	hp (W)	1 (746)
Maximum Speed ²	mph (km/hr)	5 (8)
Maximum Climbing Angle / Safe Climbing Angle(based on user weight of 90kg/200lb)	% (°)	14 (8)*
Range (per charge)(After the battery and mechanical moving parts fully break in)	mile (km)	25 (40)
Turning Radius	in (mm)	54 (1372)
Ground Clearance	in (mm)	3 (76)
Handle Bar	Type	T-Bar

SHOPRIDER

Toll-Free: 800.743.0772

www.shoprider.com