



Mobie™ Foldable Scooter

ACCESSORIES

Armrest

The Mobie Armrest is designed to provide arm comfort and support. It features a cushioned vinyl armrest that folds with the unit (no removal required when folding). The Mobie armrest can be installed in minutes.

Technical Specifications

- **Material:** Metal, foam cushion and vinyl
- **Dimensions:** 16.3" X 8" X 3.5"
- **Net Weight:** 4 lbs.
- **What's included:** Left and right armrest, mounting hardware



Easy Assembly



Battery Docking Station



The Mobie Battery Docking Station is a stand-alone portable charger enabling convenient charging of the Mobie lithium battery.

Technical Specifications

- **Part Number:** M-DS01-13
- **Material:** ABS Plastic
- **What's included:** Docking Station, Power supply
- **Color:** Black
- **Dimensions:** 8.3" X 3.5" X 2.8"
- **Net Weight:** 3 lbs.
- **Power Supply Cord Length:**
 - XLR cable: 53"
 - Power Cable: 61"
- **Certifications:** UL



Seat Cover Set

The Mobie Seat cover set is made of durable vinyl. It features a Velcro locking system and is design for easy installation.

Technical Specifications

- **Material:** Durable Vinyl
- **Color:** Black
- **What's included:** Seat backrest cover, seat cover

