

zip'r mobility
 mobility solutions corporation

zip'r3



Simplicity in Mobility

zip'r4



MODEL

ZIP'R3

ZIP'R4

BATTERY (2) 12V, 12AH DEEP CYCLE, SEALED
 CERTIFICATION CE
 CHARGER OFF BOARD, 2A
 CLEARANCE 2.1" INCH
 COLORS RED,BLUE
 DIMENSIONS 36.4" X 19.5" INCH
 HEAVIEST PIECE 29 LBS OR 13 KG
 MAX SPEED 4.25 MPH OR 6.8 KM/H
 MOTOR 155W, .21HP
 RANGE PER CHARGE 10 MILES OR 16KM
 SHIPPING WEIGHT 107 LBS OR 48 KG
 STANDARD SEAT FOLDABLE FOAM, CLOTH/VINYL
 TIRES FRONT 2.0" X 7.5" INCH
 TIRES REAR 2.5" X 7.5" INCH
 TOTAL WEIGHT 102 LBS OR 46 KG
 TURNING RADIUS 30" INCH
 WEIGHT CAPACITY 250 LBS OR 113 KG

BATTERY (2) 12V, 12AH DEEP CYCLE, SEALED
 CERTIFICATION CE
 CHARGER OFF BOARD, 2A
 CLEARANCE 2.1" INCH
 COLORS RED,BLUE
 DIMENSIONS 36.4" X 19.5" INCH
 HEAVIEST PIECE 29 LBS OR 13 KG
 MAX SPEED 4.25 MPH OR 6.8 KM/H
 MOTOR 155W, .21HP
 RANGE PER CHARGE 10 MILES OR 16KM
 SHIPPING WEIGHT 110 LBS OR 50 KG
 STANDARD SEAT FOLDABLE FOAM, CLOTH/VINYL
 TIRES FRONT 2.0" X 7.5" INCH
 TIRES REAR 2.5" X 7.5" INCH
 TOTAL WEIGHT 105 LBS OR 48 KG
 TURNING RADIUS 37" INCH
 WEIGHT CAPACITY 250 LBS OR 113 KG

To order call (425) 837-9884 or Fax (425) 837-9886 or www.zipr.com