

PINNACLE SL600HD

HEAVY DUTY STAIRLIFT



THE INDUSTRY'S ONLY 600 LB CAPACITY STAIRLIFT

High-backed, contoured extra wide seat ensures a comfortable and secure ride. Sturdy, stationary armrests and reinforced footrest add stability for safer transfers. The seat slides out from the wall as it pivots at the top and bottom landings.

As with all Pinnacle stairlifts, the HD model is engineered to the highest industry standards and designed for optimal comfort, convenience, and aesthetics.

The patented nylon polymer worm gear provides low friction for longer battery operation during a power outage and the elimination of maintenance related to track lubrication providing extra peace of mind.



UNIQUE SLIDE-SWIVEL SEAT

High capacity seat pivots 85° as it slides away from the wall to allow for installation closer to the wall and easier dismount at the top and the bottom of the stairwell.

COMFORT AND PERFORMANCE



I can safely access the upper level of my home with my stairlift. It allows me to move comfortably in and out of my chair with its swivel seat and sturdy armrests.

SL600HD PRODUCT SPECIFICATIONS

WEIGHT CAPACITY	600 lb	FOLDED WIDTH	16.5" from wall
TRACK LENGTH	15' 6" to 40'	SPEED	14 FPM
SEAT SIZE	18.5" x 14"	BACKREST SIZE	17.6" x 17"
INCLINE LIMITS	27° to 48°* * inclines above 45° may require local jurisdiction approval.	DRIVE	- Patented nylon polymer gear rack - Helical worm gear
SAFETY FEATURES	- 85° seat swivel interlock - Centrifugal overspeed brake - Footrest obstruction sensors - Lockable battery + breaker switch - Motor brake - Obstruction sensors - Seat belt	CERTIFIED	- ETL listed - CSA B44.1/ASME A17.5 - Safety standard ASME A18.1 - Rated commercial and residential
OPTIONS	- Key lock - Extra rail - Manual lowering tool	WARRANTY	- 3-year parts - 10-year gear rack - 1-year battery