



R82[®] 



rabbit[™]



Mobile standing frame for active children
mobility



The new Rabbit is suitable for children between the age of 1 to 18 years. The functionality of the Rabbit stays the same. It still provides children with the opportunity to explore and enjoy new and exciting adventures. The standing position entices the child to actively participate in things around him or her. The mobility of the Rabbit is similar to that experienced by active wheelchair users.




Chest support




Extension




Foot plate with heel stops




Knee support, hinged





Size 1



Size 2+3



Size 4


rabbit™ by R82®
Measurements

	Rabbit size 1 w/22" wheels cm (inch)	Rabbit size 2 w/26" wheels cm (inch)	Rabbit size 3 w/32" wheels cm (inch)	Rabbit size 4 w/36" wheels cm (inch)
Length, frame (L)	62 (24 1/4")	75 (29 1/4")	75 (29 1/4")	85 (33 1/4")
Length, total	62 (24 1/4")	75 (29 1/4")	80 (31 1/4")	90 (35")
Width (W)	61 1/2 (24")	68 (26 1/2")	69 (27")	70 (27 1/4")
Height from floor. Frame in highest position (H1)	100 (39")	118 (46")	131 (51")	153 (59 3/4")
Height from floor. Frame in lowest position (H1)	67 (26 1/4")	85 (33 1/4")	99 (38 1/2")	120 (46 3/4")
Height from foot plate. Frame in highest position (H2)	90 (35")	108 (42")	119 (46 1/2")	139 (54 1/4")
Height from foot plate. Frame in lowest position (H2)	57 (22 1/4")	75 (29 1/4")	86 (33 1/2")	106 (41 1/4")
Angle of the centre bar (A1)	0 to +30°	0 to +30°	0 to +30°	0 to +30°
Angle of the foot plate (A2)	+15° to -10°	+15° to -10°	+15° to -10°	+15° to -10°
Camber	8°	4°	4°	4°
Weight	13,5 kg (29,5 lb)	15,5 kg (34 lb)	16 kg (35 lb)	18,5 kg (40,5 lb)
Max. load/user weight	40 kg (88 lb)	50 kg (110 lb)	60 kg (132 lb)	70 kg (154 lb)

R82 A/S

Parallelsvej 3
8751 Gedved
Denmark
Tel +45 79 68 58 88
Fax +45 75 66 51 92
E-mail: r82@r82.com
Web: www.r82.com

R82 UK Ltd.

Unit D4A
Coombswood Business Park East
Coombswood Way, Halesowen
West Midlands B62 8BH
Great Britain
Tel 0121 561 2222
Fax 0121 559 5437
E-mail: info@r82.com
Web: www.r82.com

Distributor: